



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Esophagus (normal) FFPE tissue slide

Catalog Number: S1307

Regulatory Status: For research use only (RUO)

Product Description: Esophagus (normal) FFPE tissue slide contains tissue section from a single Rhesus (Macaca Mulatta) donor.

Tissue: Esophagus

Host: Macaque

Clinical Diagnosis: Normal

Recommended Reagents: [IHC antibodies](#)

[CISH probes](#)

[FISH probes](#)

[mutaFISH™ probes](#)

Videos: [H&E Staining](#)

[Immunohistochemistry](#)

Services: [FFPE Tumor Tissue DNA Service](#)

[Immunohistochemistry Service](#)

[Custom FISH Probe Service](#)

[Custom mutaFISH™ Probe Service](#)

Recommend Reagents: [IHC antibodies](#)

[CISH probes](#)

[FISH probes](#)

[mutaFISH™ probes](#)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Tissue sections are 5 um in thickness, and mounted on positively charged glass slides.

Recommend Usage: Please bake slides at 60°C for 30 minutes before use.

Storage Instruction: Store at 4°C.