



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

GST Selector

Catalog Number: U0390

Regulation Status: For research use only (RUO)

Product Description: GST Selector is based on a high-affinity single-domain antibody (sdAb) that is compatible not only with physiological buffers but also with high stringency buffers.

Clone Name: G19

Host: Camelid

Applications: ChIP, Co-IP, IP
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Recognizes *Schistosoma japonicum* Glutathione S-transferase (GST).

Form: Liquid

Recommend Usage: ChIP
Co-Immunoprecipitation
Immunoprecipitation
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (20% ethanol, 50% slurry).

Storage Instruction: Store at 4°C. Do not freeze.