



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### GFP Selector

**Catalog Number:** U0391

**Regulation Status:** For research use only (RUO)

**Product Description:** GFP Selector is based on a high-affinity single-domain antibody (sdAb) that is compatible not only with physiological buffers but also with high stringency buffers.

**Clone Name:** 1H1

**Host:** Camelid

**Applications:** ChIP, Co-IP, IP  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Recognizes GFP, mEGFP, superfolder GFP, Emerald GFP and most common CFP and YFP variants. It does not cross-react with mCherry, mRFP, dsRed, mTagBFP or their most common derivatives.

**Form:** Liquid

**Recommend Usage:** ChIP  
Co-Immunoprecipitation  
Immunoprecipitation  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (20% ethanol, 50% slurry).

**Storage Instruction:** Store at 4°C. Do not freeze.