



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Forensic DNA Purification Kit (Pre-filled)

Catalog Number: U0432

Regulation Status: For research use only (RUO)

Product Description: Forensic DNA Purification Kit (Pre-filled) is designed for the optimal automated extraction of DNA from a broad range of forensic samples, including human blood stain, dried blood spot, hair follicle, semen, cigarette, chewing gum, and nail.

Pack Size: 96 reactions x 4

Suitable Sample: Forensic samples (blood stain, dried blood spot, hair follicle, semen, cigarette, chewing gum, nail)

Supplied Product: Kit content:

1. Auto Plate
 - a. Lysis Buffer
 - a. Wash Buffer 1
 - b. Wash Buffer 2
 - c. Magnetic Beads
 - d. Elution Buffer
2. Proteinase K
3. Incubation Buffer
4. Elution Buffer
5. Spin tips

Instrument: [Spiralator™ 48](#)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Instruction: Components under room temperature (15~35°C) can be stored until the expiration date labeled on the box. The Proteinase K is transported at room temperature. Upon received, please store Proteinase K at 4°C.