



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### Stool DNA Purification Kit (Pre-filled)

**Catalog Number:** U0433

**Regulation Status:** For research use only (RUO)

**Product Description:** Stool DNA Purification Kit (Pre-filled) is designed for the optimal automated extraction of DNA from stool of human and other species.

**Pack Size:** 96 reactions x 4

**Suitable Sample:** Stool

**Sample Size:** 50 mg stool

**Typical DNA yield:** 2 ug

**Supplied Product:** Kit content:

1. Auto Plate
  - a. Lysis Buffer
  - a. Wash Buffer 1
  - b. Wash Buffer 2
  - a. Wash Buffer 3
  - c. Magnetic Beads
  - d. Elution Buffer
2. Proteinase K
3. Elution Buffer
4. Incubation Buffer-B845
5. Incubation Buffer-B871
6. Spin tips

**Instrument:** [Spiralator™ 48](#)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Instruction:** Components under room temperature (15~35°C) can be stored until the expiration date labeled on the box. The Proteinase K is transported at room temperature. Upon received, please store Proteinase K at 4°C.