



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BC05023a : Colon cancer tissue array with matched cancer adjacent tissue, including TNM and pathology grade, 18 cases/ 54 cores, replacing BC05023

Microarray Panel	Colon carcinoma tissue microarray with matched cancer adjacent tissue, containing 18 cases of adenocarcinoma, triplicate cores per case		
Cores	54	BC05023a H&E 	
Cases	18		
Row number	6		
Column number	9		
Core Diameter (mm)	1.5		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.		

	1	2	3	4	5	6	7	8	9
A	Col	Col	Col	Col	Col	Col	Col	Col	Col
B	Col	Col	Col	Col	Col	Col	Col	Col	Col
C	Col	Col	Col	Col	Col	Col	Col	Col	Col
D	Col	Col	Col	Col	Col	Col	Col	Col	Col
E	Col	Col	Col	Col	Col	Col	Col	Col	Col
F	Col	Col	Col	Col	Col	Col	Col	Col	Ad*

Legend: Col - Colon
 - AT, - Malignant tumor, - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	32	F	Colon	Adenocarcinoma	T3N0M0	1		Malignant	
A2	2	32	F	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-		AT	
A3	3	32	F	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-		NAT	
A4	4	42	M	Colon	Adenocarcinoma	T3N1M0	1		Malignant	
A5	5	42	M	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-		AT	
A6	6	42	M	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-		NAT	
A7	7	55	F	Colon	Mucinous adenocarcinoma	T3N1M0	1		Malignant	
A8	8	55	F	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-		AT	
A9	9	55	F	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-		NAT	
B1	10	50	M	Colon	Adenocarcinoma	T3N0M0	1		Malignant	
B2	11	50	M	Colon	Cancer adjacent colon tissue	-	-		AT	
B3	12	50	M	Colon	Adjacent normal colon tissue	-	-		NAT	
B4	13	68	F	Colon	Adenocarcinoma	T4N2M0	1		Malignant	
B5	14	68	F	Colon	Cancer adjacent colon tissue	-	-		AT	
B6	15	68	F	Colon	Adjacent normal colon tissue	-	-		NAT	
B7	16	62	M	Colon	Adenocarcinoma	T4N0M0	1--2		Malignant	
B8	17	62	M	Colon	Cancer adjacent colon tissue	-	-		AT	
B9	18	62	M	Colon	Adjacent normal colon tissue	-	-		NAT	
C1	19	62	F	Colon	Adenocarcinoma	T4N1M0	1--2		Malignant	
C2	20	62	F	Colon	Cancer adjacent colon tissue	-	-		AT	
C3	21	62	F	Colon	Adjacent normal colon tissue	-	-		NAT	
C4	22	59	F	Colon	Adenocarcinoma	T3N1M0	2		Malignant	
C5	23	59	F	Colon	Cancer adjacent colon tissue	-	-		AT	

C6	24	59	F	Colon	Adjacent normal colon tissue	-	-	NAT	
C7	25	46	F	Colon	Mucinous adenocarcinoma	T3N0M0	2	Malignant	
C8	26	46	F	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-	AT	
C9	27	46	F	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-	NAT	
D1	28	52	F	Colon	Adenocarcinoma	T3N1M0	2	Malignant	
D2	29	52	F	Colon	Cancer adjacent colon tissue	-	-	AT	
D3	30	52	F	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-	NAT	
D4	31	33	F	Colon	Adenocarcinoma	T4N1M0	2	Malignant	
D5	32	33	F	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-	AT	
D6	33	33	F	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-	NAT	
D7	34	61	M	Colon	Adenocarcinoma	T4N0M0	1	Malignant	
D8	35	61	M	Colon	Cancer adjacent colon tissue	-	-	AT	
D9	36	61	M	Colon	Adjacent normal colon tissue	-	-	NAT	
E1	37	67	M	Colon	Adenocarcinoma	T3N1M0	2	Malignant	
E2	38	67	M	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-	AT	
E3	39	67	M	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-	NAT	
E4	40	53	F	Colon	Adenocarcinoma	T3N0M0	-	Malignant	
E5	41	53	F	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-	AT	
E6	42	53	F	Colon	Adjacent normal colon tissue	-	-	NAT	
E7	43	53	F	Colon	Adenocarcinoma	T3N0M0	3	Malignant	
E8	44	53	F	Colon	Cancer adjacent colon tissue	-	-	AT	
E9	45	53	F	Colon	Adjacent normal colon tissue	-	-	NAT	
F1	46	34	F	Colon	Adenocarcinoma	T2N0M0	2	Malignant	
F2	47	34	F	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-	AT	
F3	48	34	F	Colon	Adjacent normal colon tissue	-	-	NAT	
F4	49	74	F	Colon	Adenocarcinoma	T3N1M0	3	Malignant	
F5	50	74	F	Colon	Cancer adjacent colon tissue	-	-	AT	
F6	51	74	F	Colon	Adjacent normal colon tissue	-	-	NAT	
F7	52	52	F	Colon	Adenocarcinoma	T3N1M0	3	Malignant	
F8	53	52	F	Colon	Cancer adjacent colon tissue (chronic inflammation of mucous membrane)	-	-	AT	
F9	54	52	F	Colon	Adjacent normal colon tissue (chronic inflammation of mucous membrane)	-	-	NAT	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)	-	-	Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada



1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU



Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World



184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland



Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland



T: +49 (0) 69 779099
F: +49 (0) 69 13376880