



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BC09022 : Endometrial carcinoma tissue array with cancer adjacent tissue and adjacent normal tissue, including TNM, clinical stage and pathology grade, 67 cases/72 cores

Microarray Panel	Endometrial carcinoma tissue microarray, containing 24 each cases of adjacent normal tissue, endometrial carcinoma and cancer adjacent tissue, single core per case	
Cores	72	
Cases	67	
Row number	8	
Column number	9	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “**” means invalid core; “-” means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8	9
A	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute
B	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute
C	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute
D	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute
E	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute
F	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute
G	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute
H	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute	Ute

Legend: Ute - Uterus

● - AT, ● - Malignant tumor, ● - Malignant tumor (stage I), ● - Malignant tumor (stage IB), ● - Malignant tumor (stage IC), ● - Malignant tumor (stage II), ● - Malignant tumor (stage IIB), ● - Malignant tumor (stage IIIC), ● - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	62	F	Uterus	Senile endometritis	-	-	-	NAT	
A2	2	59	F	Uterus	Endometrioid adenocarcinoma	T1aN0M0	2	I	Malignant	
A3	3	48	F	Uterus	Atrophic endometrium	-	-	-	AT	
A4	4	57	F	Uterus	Adjacent normal endometrium tissue	-	-	-	NAT	
A5	5	62	F	Uterus	Endometrioid adenocarcinoma	T1N0M0	2	I	Malignant	
A6	6	43	F	Uterus	Cancer adjacent endometrium tissue	-	-	-	AT	
A7	7	42	F	Uterus	Endometritis	-	-	-	NAT	
A8	8	42	F	Uterus	Endometrioid adenocarcinoma (uterine smoot	-	-	-	Malignant	

G6	60	55	F	Uterus	Endometritis with atypical hyperplasia of epithelium	-	-	-	AT	
G7	61	60	F	Uterus	Adjacent normal endometrium tissue (uterine smooth muscle tissue)	-	-	-	NAT	
G8	62	60	F	Uterus	Endometrioid adenocarcinoma	T1cN0M0	3	IC	Malignant	
G9	63	56	F	Uterus	Cancer adjacent endometrium tissue	-	-	-	AT	
H1	64	69	F	Uterus	Senile endometrium	-	-	-	NAT	
H2	65	69	F	Uterus	Endometrioid adenocarcinoma	T1N0M0	3	I	Malignant	
H3	66	52	F	Uterus	Cancer adjacent endometrium tissue	-	-	-	AT	
H4	67	57	F	Uterus	Adjacent normal endometrium tissue	-	-	-	NAT	
H5	68	57	F	Uterus	Endometrioid adenocarcinoma	T1cN0M0	3	IC	Malignant	
H6	69	32	F	Uterus	Cancer adjacent endometrium tissue	-	-	-	AT	
H7	70	60	F	Uterus	Adjacent normal endometrium tissue (uterine smooth muscle tissue)	-	-	-	NAT	
H8	71	64	F	Uterus	Endometrioid adenocarcinoma (endometritis with atypical hyperplasia of epithelium)	T1cN0M0		IC	Malignant	
H9	72	57	F	Uterus	Cancer adjacent endometrium tissue (uterine smooth muscle tissue)	-	-	-	AT	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

TNM grading:

T - Primary tumor

- Tx - Primary tumor cannot be assessed
- T0 - No evidence of primary tumor
- Tis - Carcinoma in situ; intraepithelial or invasion of lamina propria
- T1 - Tumor invades submucosa
- T2 - Tumor invades muscularis propria
- T3 - Tumor invades through muscularis propria into subserosa or into non-peritonealized pericolic or perirectal tissues.
- T4 - Tumor directly invades other organs or structures and/or perforate visceral peritoneum

N - Regional lymph nodes

- Nx - Regional lymph nodes cannot be assessed
- N0 - No regional lymph node metastasis
- N1 - Metastasis in 1 to 3 regional lymph nodes
- N2 - Metastasis in 4 or more regional lymph nodes

M - Distant metastasis

- Mx - Distant metastasis cannot be assessed
- M0 - No distant metastasis
- M1 - Distant metastasis

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, U.K.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Via Lisano 3
(CP.683)
6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85