



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



BO244g : Normal bone tissue array with osteosarcoma tissue, including TNM and clinical stage, 12 cases/24 cores, replacing BO244f

Microarray Panel	Normal bone tissue microarray with osteosarcoma tissue, containing 11 cases of bone tissue, 1 osteosarcoma, duplicate cores per case	
Cores	24	
Cases	12	
Row number	3	
Column number	8	
Core Diameter (mm)	2	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8
A	Bon	Bon	Bon	Bon	Bon	Bon	Bon	Bon
B	Bon	Bon	Bon	Bon	Bon	Bon	Bon	Bon
C	Bon	Bon	Bon	Bon	Bon	Bon	Bon	Adn

Legend: Bon - Bone/Left tibia
● - Malignant tumor, ● - Malignant tumor (stage IIB), ● - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	58	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
A2	2	58	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
A3	3	57	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
A4	4	57	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
A5	5	56	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
A6	6	56	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
A7	7	56	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
A8	8	56	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B1	9	62	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B2	10	62	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B3	11	70	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B4	12	70	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B5	13	75	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B6	14	75	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B7	15	67	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
B8	16	67	M	Bone/Rib	Bone tissue	-		-	NAT	
C1	17	58	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
C2	18	58	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	
C3	19	41	M	Bone/Rib	Bone and bone marrow tissue	-		-	NAT	

C4	20	41	M	Bone/Rib	Bone tissue	-	-	NAT	
C5	21	66	F	Bone/Rib	Bone and bone marrow tissue	-	-	NAT	
C6	22	66	F	Bone/Rib	Bone and bone marrow tissue	-	-	NAT	
C7	23	19	M	Bone/Left tibia	Osteosarcoma of left tibia	T2N0M0 G3	IIB	Malignant	
C8	24	19	M	Bone/Left tibia	Osteosarcoma of left tibia	T2N0M0 G3	IIB	Malignant	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-	Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
 184 Park Drive, Milton Park
 Abingdon OX14 4SE.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

North America 
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

Europe 
 Berenkoog 41,
 1822 BH Alkmaar,
 Netherlands
 T: +31 (0) 72 8080244
 F: +31 (0) 72 8080142

Switzerland 
 Via Lisano 3,
 (CP.683)
 CH-6900
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Germany 
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880