



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BR2410 : Breast invasive lobular carcinoma tissue array, including TNM, clinical stage and pathology grade, 24 cases/ 24 cores

Microarray Panel	Breast invasive lobular carcinoma tissue microarray, containing 24 cases of invasive lobular carcinoma, single core per case		
Cores	24		
Cases	24		
Row number	4		
Column number	6		
Core Diameter (mm)	1.5		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “*” means invalid core; “-” means no applicable or negative in IHC markers.		

	1	2	3	4	5	6
A	Bre	Bre	Bre	Bre	Bre	Bre
B	Bre	Bre	Bre	Bre	Bre	Bre
C	Bre	Bre	Bre	Bre	Bre	Bre
D	Bre	Bre	Bre	Bre	Bre	Adr

Legend: Bre - Breast

● - Malignant tumor, ● - Malignant tumor (stage IIA), ● - Malignant tumor (stage IIB), ● - Malignant tumor (stage IIIA), ● - Malignant tumor (stage IIIB)

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	59	F	Breast	Invasive lobular carcinoma (cataplasia tissue)	T2N0M0		IIA	malignant	
A2	2	47	F	Breast	Invasive lobular carcinoma	T2N1M0		IIB	malignant	
A3	3	37	F	Breast	Invasive lobular carcinoma	T4N1M0		IIIB	malignant	
A4	4	49	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
A5	5	73	F	Breast	Invasive lobular carcinoma	T4N0M0		IIIB	malignant	
A6	6	35	F	Breast	Invasive lobular carcinoma	T3N0M0		IIB	malignant	
B1	7	46	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
B2	8	51	F	Breast	Invasive lobular carcinoma	T4N0M0		IIIB	malignant	
B3	9	19	F	Breast	Invasive lobular carcinoma	T2N1M0		IIB	malignant	
B4	10	37	F	Breast	Invasive lobular carcinoma (cataplasia tissue)	T2N0M0		IIA	malignant	
B5	11	45	F	Breast	Invasive lobular carcinoma	T2N1M0		IIB	malignant	
B6	12	50	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
C1	13	41	F	Breast	Invasive lobular carcinoma	T4N0M0		IIIB	malignant	
C2	14	53	F	Breast	Invasive lobular carcinoma (sparse)	T4N1M0		IIIB	malignant	
C3	15	47	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
C4	16	49	F	Breast	Invasive lobular carcinoma	T2N1M0		IIB	malignant	
C5	17	69	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
C6	18	58	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
D1	19	63	F	Breast	Invasive lobular carcinoma	T3N0M0		IIB	malignant	
D2	20	50	F	Breast	Invasive lobular carcinoma	T2N2M0		IIIA	malignant	

D3	21	42	F	Breast	Invasive lobular carcinoma	T2N1M0		IIB	malignant	
D4	22	60	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
D5	23	57	F	Breast	Invasive lobular carcinoma	T2N0M0		IIA	malignant	
D6	24	50	F	Breast	Invasive lobular carcinoma	T2N1M0		IIB	malignant	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada



1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU



Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World



184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland



Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland



T: +49 (0) 69 779099
F: +49 (0) 69 13376880