



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BR246d : Breast cancer with breast tissue array, including pathology grade, TNM and clinical stage (reference AJCC 7th version), 12 cases/24 cores, replacing BR246c

Microarray Panel	Breast invasive carcinoma with breast tissue microarray, containing 10 cases of invasive carcinoma of no special type, 1 each of cancer adjacent breast tissue and adjacent normal breast tissue, duplicated cores per case		
Cores	24	BR246d H&E 	
Cases	12		
Row number	4		
Column number	6		
Core Diameter (mm)	1.5		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.		

	1	2	3	4	5	6
A	Bre	Bre	Bre	Bre	Bre	Bre
B	Bre	Bre	Bre	Bre	Bre	Bre
C	Bre	Bre	Bre	Bre	Bre	Bre
D	Bre	Bre	Bre	Bre	Bre	Adr

Legend: Bre - Breast

● - AT, ● - Malignant tumor, ● - Malignant tumor (stage IIA), ● - Malignant tumor (stage IIB), ● - Malignant tumor (stage IIIA), ● - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	45	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
A2	2	45	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
A3	3	42	F	Breast	Invasive carcinoma of no special type	T2N2M0	2	IIIA	Malignant	
A4	4	42	F	Breast	Invasive carcinoma of no special type	T2N2M0	2	IIIA	Malignant	
A5	5	48	F	Breast	Invasive carcinoma of no special type	T2M0N0	2	IIA	Malignant	
A6	6	48	F	Breast	Invasive carcinoma of no special type	T2M0N0	2	IIA	Malignant	
B1	7	56	F	Breast	Invasive carcinoma of no special type	T2N1M0	2	IIB	Malignant	
B2	8	56	F	Breast	Invasive carcinoma of no special type	T2N1M0	2	IIB	Malignant	
B3	9	44	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
B4	10	44	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
B5	11	47	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
B6	12	47	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
C1	13	47	F	Breast	Invasive carcinoma of no special type (fibrous tissue and fatty tissue)	T2N0M0	*	IIA	Malignant	
C2	14	47	F	Breast	Invasive carcinoma of no special type	T2N0M0	3	IIA	Malignant	
C3	15	48	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
C4	16	48	F	Breast	Invasive carcinoma of no special type	T2N0M0	2	IIA	Malignant	
C5	17	29	F	Breast	Invasive carcinoma of no special type	T2N0M0	3	IIA	Malignant	
C6	18	29	F	Breast	Invasive carcinoma of no special type	T2N0M0	3	IIA	Malignant	
D1	19	52	F	Breast	Invasive carcinoma of no special type	T2N0M0	3	IIA	Malignant	
D2	20	52	F	Breast	Invasive carcinoma of no special type	T2N0M0	3	IIA	Malignant	
D3	21	45	F	Breast	Adjacent normal breast tissue	-	-	-	NAT	
D4	22	45	F	Breast	Adjacent normal breast tissue	-	-	-	NAT	
D5	23	54	F	Breast	Cancer adjacent breast tissue	-	-	-	AT	
D6	24	54	F	Breast	Cancer adjacent breast tissue	-	-	-	AT	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)	-	-	-	Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland
Via Lisano 3,
6900 Masnago (CP 683)
CH-6900 Masnago
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland
T: +49 (0) 69 779099
F: +49 (0) 69 13376880