



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

CO241a : Colon cancer tissue array with adjacent normal colon tissue, including TNM, clinical stage and pathology grade, 6 cases/24 cores

Microarray Panel	Colonic carcinoma tissue microarray with adjacent normal colon tissue, containing 6 colon adenocarcinoma, quadruple cores per case	
Cores	24	CO241a H&E 
Cases	6	
Row number	3	
Column number	8	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Quality Control	Anti-Actin confirmed	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “*” means invalid core; “-” means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8
A	Col	Col	Col	Col	Col	Col	Col	Col
B	Col	Col	Col	Col	Col	Col	Col	Col
C	Col	Col	Col	Col	Col	Col	Col	Col

Legend: Col - Colon

● - Malignant tumor, ● - Malignant tumor (stage IIIc), ● - Malignant tumor (stage IIa), ● - Malignant tumor (stage IIb), ● - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	62	F	Colon	Adenocarcinoma	T4N0M0	1	IIb	malignant	
A2	2	62	F	Colon	Adenocarcinoma	T4N0M0	1	IIb	malignant	
A3	3	62	F	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
A4	4	62	F	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
A5	5	39	M	Colon	Adenocarcinoma	T4N0M0	1	IIb	malignant	
A6	6	39	M	Colon	Adenocarcinoma	T4N0M0	1	IIb	malignant	
A7	7	39	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
A8	8	39	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
B1	9	62	M	Colon	Adenocarcinoma	T3N0M0	2	IIa	malignant	
B2	10	62	M	Colon	Adenocarcinoma	T3N0M0	2	IIa	malignant	
B3	11	62	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
B4	12	62	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
B5	13	42	F	Colon	Adenocarcinoma	T3N0M0	2	IIa	malignant	
B6	14	42	F	Colon	Adenocarcinoma	T3N0M0	2	IIa	malignant	
B7	15	42	F	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
B8	16	42	F	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
C1	17	61	M	Colon	Adenocarcinoma	T3N2M0	2	IIIc	malignant	
C2	18	61	M	Colon	Adenocarcinoma	T3N2M0	2	IIIc	malignant	

C3	19	61	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
C4	20	61	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
C5	21	56	M	Colon	Adenocarcinoma	T4N2M0	3	IIIc	malignant	
C6	22	56	M	Colon	Adenocarcinoma	T4N2M0	3	IIIc	malignant	
C7	23	56	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
C8	24	56	M	Colon	Cancer adjacent normal colon tissue	-	-	-	NAT	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880