



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

CO802d : Colon cancer tissue array with adjacent normal tissue, including pathology grade, TNM and clinical stage, 80 cases/80 cores, replacing CO802c

Microarray Panel	Colon adenocarcinoma tissue microarray with adjacent normal tissue, containing 71 cases of adenocarcinoma, 5 signet-ring cell carcinoma, 2 squamous cell carcinoma, plus 2 adjacent normal colon tissue, single core per case	
Cores	80	
Cases	80	
Row number	8	
Column number	10	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8	9	10
A	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col
B	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col
C	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col
D	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col
E	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col
F	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col
G	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col
H	Col	Col	Col	Col	Col	Col	Col	Col	Col	Col

Legend: Col - Colon

● - Malignant tumor, ● - Malignant tumor (stage I), ● - Malignant tumor (stage IIA), ● - Malignant tumor (stage IIB), ● - Malignant tumor (stage IIIB), ● - Malignant tumor (stage IIIC), ● - Malignant tumor (stage IV), ● - NAT, ● - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	63	F	Colon	Tubular adenocarcinoma	T2N0M0	1	I	Malignant	
A2	2	53	M	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
A3	3	55	M	Colon	Tubular adenocarcinoma	T3N2M0	1	IIIB	Malignant	
A4	4	61	M	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
A5	5	48	M	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
A6	6	54	M	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
A7	7	49	F	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
A8	8	63	F	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
A9	9	31	M	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	
A10	10	32	M	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	
B1	11	66	M	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	
B2	12	62	M	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	
B3	13	61	M	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	

B4	14	62	F	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	
B5	15	52	F	Colon	Tubular adenocarcinoma	T3N1M0	1	IIIB	Malignant	
B6	16	38	M	Colon	Tubular adenocarcinoma	T3N1M0	2	IIIB	Malignant	
B7	17	61	M	Colon	Tubular adenocarcinoma	T3N1M0	1	IIIB	Malignant	
B8	18	62	F	Colon	Tubular adenocarcinoma	T4N1M0	1	IIIB	Malignant	
B9	19	54	F	Colon	Tubular adenocarcinoma	T3N0M0	2	IIA	Malignant	
B10	20	50	F	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
C1	21	45	F	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
C2	22	46	M	Colon	Tubular adenocarcinoma	T3N0M0	-	IIA	Malignant	
C3	23	64	M	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	
C4	24	38	M	Colon	Tubular adenocarcinoma	T4N1M0	1	IIIB	Malignant	
C5	25	51	M	Colon	Tubular adenocarcinoma	T4N0M0	1	IIB	Malignant	
C6	26	63	M	Colon	Tubular adenocarcinoma	T2N0M0	1	I	Malignant	
C7	27	51	F	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
C8	28	71	M	Colon	Tubular adenocarcinoma	T3N0M0	2	IIA	Malignant	
C9	29	73	M	Colon	Papillary tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
C10	30	55	F	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
D1	31	54	M	Colon	Tubular adenocarcinoma	T3N0M0	2	IIA	Malignant	
D2	32	52	M	Colon	Tubular adenocarcinoma	T3N0M0	2	IIA	Malignant	
D3	33	76	M	Colon	Tubular adenocarcinoma	T3N0M0	2	IIA	Malignant	
D4	34	53	F	Colon	Tubular adenocarcinoma	T3N0M0	2	IIA	Malignant	
D5	35	49	M	Colon	Tubular adenocarcinoma	T3N0M0	1	IIA	Malignant	
D6	36	66	M	Colon	Tubular adenocarcinoma	T4N0M0	2	IIB	Malignant	
D7	37	73	M	Colon	Tubular adenocarcinoma	T4N0M0	2	IIB	Malignant	
D8	38	61	F	Colon	Tubular adenocarcinoma	T3N1M0	2	IIIB	Malignant	
D9	39	58	F	Colon	Tubular adenocarcinoma	T3N1M0	2	IIIB	Malignant	
D10	40	50	M	Colon	Papillary adenocarcinoma	T3N0M0	1	IIA	Malignant	
E1	41	61	M	Colon	Papillary adenocarcinoma	T2N0M0	1	I	Malignant	
E2	42	60	M	Colon	Papillary adenocarcinoma	T2N0M0	1	I	Malignant	
E3	43	49	F	Colon	Papillary adenocarcinoma	T3N0M0	1	IIA	Malignant	
E4	44	71	M	Colon	Papillary adenocarcinoma	T3N0M0	1	IIA	Malignant	
E5	45	32	M	Colon	Tubular adenocarcinoma	T3N1M0	1	IIIB	Malignant	
E6	46	41	M	Colon	Papillary adenocarcinoma	T4N0M0	1	IIB	Malignant	
E7	47	12	F	Colon	Papillary adenocarcinoma	T3N0M0	1	IIA	Malignant	
E8	48	50	F	Colon	Tubular adenocarcinoma	T4N0M0	-	IIB	Malignant	
E9	49	74	M	Colon	Tubular adenocarcinoma	T3N1M0	2	IIIB	Malignant	
E10	50	71	M	Colon	Tubular adenocarcinoma	T3N0M0	2	IIA	Malignant	
F1	51	67	F	Colon	Adenocarcinoma	T3N0M0	2	IIA	Malignant	
F2	52	56	F	Colon	Adenocarcinoma	T3N0M0	3	IIA	Malignant	
F3	53	72	M	Colon	Mucinous adenocarcinoma	T3N0M0	3	IIA	Malignant	
F4	54	65	M	Colon	Adenocarcinoma	T4N0M0	3	IIB	Malignant	
F5	55	20	F	Colon	Adenocarcinoma	T4N1M0	3	IIIB	Malignant	
F6	56	48	M	Colon	Adenocarcinoma	T4N2M0	3	IIIC	Malignant	
F7	57	49	F	Colon	Mucinous adenocarcinoma	T2N0M0	1	I	Malignant	
F8	58	68	M	Colon	Mucinous adenocarcinoma	T2N0M0	1	I	Malignant	
F9	59	51	M	Colon	Mucinous adenocarcinoma	T2N0M0	1	I	Malignant	
F10	60	44	F	Colon	Mucinous adenocarcinoma	T3N0M0	1	IIA	Malignant	
G1	61	39	M	Colon	Mucinous adenocarcinoma	T3N0M0	1	IIA	Malignant	
G2	62	70	F	Colon	Mucinous adenocarcinoma	T3N0M0	1	IIA	Malignant	
G3	63	41	F	Colon	Mucinous adenocarcinoma	T3N0M0	1	IIA	Malignant	
G4	64	50	M	Colon	Mucinous adenocarcinoma	T3N0M0	1	IIA	Malignant	
G5	65	45	F	Colon	Mucinous adenocarcinoma	T4N2M0	1	IIIC	Malignant	
G6	66	68	F	Colon	Mucinous adenocarcinoma	T3N1M0	1	IIIB	Malignant	
G7	67	74	F	Colon	Mucinous adenocarcinoma	T3N0M0	2	IIA	Malignant	
G8	68	54	F	Colon	Mucinous adenocarcinoma	T3N0M0	2	IIA	Malignant	
G9	69	62	M	Colon	Mucinous adenocarcinoma	T3N0M0	2	IIA	Malignant	
G10	70	39	M	Colon	Mucinous adenocarcinoma	T3N1M0	3	IIIB	Malignant	
H1	71	57	M	Colon	Mucinous adenocarcinoma	T3N1M0	-	IIIB	Malignant	
H2	72	37	F	Colon	Signet-ring cell carcinoma	T3N0M0	-	IIA	Malignant	
H3	73	51	F	Colon	Signet-ring cell carcinoma	T3N0M0	-	IIA	Malignant	

H4	74	46	M	Colon	Signet-ring cell carcinoma	T3N0M0	-	IIA	Malignant	
H5	75	43	M	Colon	Signet-ring cell carcinoma	T3N0M0	-	IIA	Malignant	
H6	76	81	F	Colon	Signet-ring cell carcinoma	T4N2M0	-	IIIC	Malignant	
H7	77	56	F	Colon	Squamous cell carcinoma	T4N1M1	3	IV	Malignant	
H8	78	68	M	Colon	Squamous cell carcinoma	T4N1M0	3	IIIB	Malignant	
H9	79	57	M	Colon	Adjacent normal colon tissue	T2N0M0	-	I	NAT	
H10	80	69	F	Colon	Adjacent normal colon tissue	T4N0M0	-	IIB	NAT	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
 184 Park Drive, Milton Park
 Abingdon OX14 4SE.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

North America 
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

Europe 
 Berenkoog 41,
 1822 BH Alkmaar,
 Netherlands
 T: +31 (0) 72 8080244
 F: +31 (0) 72 8080142

Switzerland 
 Via Lisano 3,
 (CP.683)
 CH-6900
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Germany 
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880