



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet

DNA source information: Human Uterus Tumor FFPE tissue sections

Shipping Condition: Blue ice. **Storage Condition:** 2-8°C, -20°C for long term storage

Age of block: 6 years

Important

DNA is stored in Qiagen ATE Buffer.

Description

DNA is isolated using QIAamp DNA FFPE Tissue Kit (Cat. No. 56404).

Paraffin Tissue Section Name	Lot Number	DNA Qty. (µg)	Total unit(s)	Concentration (µg/µl)	A260/280	A260/230
Human Tumor: Uterus	B409140	5	1	0.2	1.95	2.22

Donor Information

Tissue Type	Pathological Diagnosis	Age	Gender
Human Tumor: Uterus	Endometrium Adenocarcinoma	58	F

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: 

UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

United States

23591 El Toro Rd, Suite #167
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703