



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

EN801b : Endocrine cancer tissue array, including TNM and clinical stage, 40 cases/80 cores, replacing EN801a

Microarray Panel	Endocrine carcinoma tissue microarray, containing 40 cases of endocrine cancer (10 adrenal cortical adenocarcinoma, 4 pancreas islet cell carcinoma, 2 each of esophagus and cardia neuroendocrine carcinoma, 3 stomach carcinoid, 1 each of colon neuroendocrine, small intestine carcinoid and rectum carcinoid, 12 lung atypical carcinoid and 4 lung carcinoid), duplicate cores per case	
Cores	80	EN801b H&E 
Cases	40	
Row number	8	
Column number	10	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8	9	10
A	Adr	Adr	Adr	Adr	Adr	Adr	Adr	Adr	Adr	Adr
B	Adr	Adr	Adr	Adr	Adr	Adr	Adr	Adr	Adr	Adr
C	Pan	Pan	Pan	Pan	Pan	Pan	Pan	Pan	Eso	Eso
D	Eso	Eso	Car	Car	Car	Car	Sto	Sto	Sto	Sto
E	Sto	Sto	Col	Col	Sma	Sma	Rec	Rec	Lun	Lun
F	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun
G	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun
H	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Lun	Adr

Legend: Adr - Adrenal gland, Car - Cardia, Col - Colon, Eso - Esophagus, Lun - Lung, Pan - Pancreas, Rec - Rectum, Sma - Small intestine, Sto - Stomach

● - Malignant tumor, ● - Malignant tumor, ● - Malignant tumor (stage I), ● - Malignant tumor (stage IB), ● - Malignant tumor (stage IIA), ● - Malignant tumor (stage IIB), ● - Malignant tumor (stage IIIB)

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	44	M	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A2	2	44	M	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A3	3	48	F	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A4	4	48	F	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A5	5	53	F	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A6	6	53	F	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A7	7	45	M	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A8	8	45	M	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A9	9	10	F	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
A10	10	10	F	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
B1	11	45	M	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	
B2	12	45	M	Adrenal gland	Adrenal cortical adenocarcinoma	-		-	Malignant	

B3	13	58	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
B4	14	58	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
B5	15	35	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
B6	16	35	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
B7	17	59	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
B8	18	59	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
B9	19	54	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
B10	20	54	M	Adrenal gland	Adrenal cortical adenocarcinoma	-	-	Malignant	
C1	21	52	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C2	22	52	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C3	23	15	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C4	24	15	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C5	25	40	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C6	26	40	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C7	27	16	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C8	28	16	F	Pancreas	Islet cell carcinoma	-	-	Malignant	
C9	29	53	F	Esophagus	Neuroendocrine carcinoma with necrosis	T2N0M0	IB	Malignant	
C10	30	53	F	Esophagus	Neuroendocrine carcinoma	T2N0M0	IB	Malignant	
D1	31	55	M	Esophagus	Neuroendocrine carcinoma	T1N0M0	IB	Malignant	
D2	32	55	M	Esophagus	Neuroendocrine carcinoma	T1N0M0	IB	Malignant	
D3	33	68	M	Cardia	Neuroendocrine carcinoma	T3N1M0	IIIB	Malignant	
D4	34	68	M	Cardia	Neuroendocrine carcinoma	T3N1M0	IIIB	Malignant	
D5	35	68	M	Cardia	Neuroendocrine carcinoma	T3N1M0	IIIB	Malignant	
D6	36	68	M	Cardia	Neuroendocrine carcinoma	T3N1M0	IIIB	Malignant	
D7	37	50	M	Stomach	Carcinoid	T4N0M0	IIB	Malignant	
D8	38	50	M	Stomach	Carcinoid	T4N0M0	IIB	Malignant	
D9	39	64	F	Stomach	Carcinoid	T2N0M0	I	Malignant	
D10	40	64	F	Stomach	Carcinoid	T2N0M0	I	Malignant	
E1	41	52	F	Stomach	Carcinoid	T2N0M0	I	Malignant	
E2	42	52	F	Stomach	Carcinoid	T2N0M0	I	Malignant	
E3	43	63	F	Colon	Neuroendocrine carcinoma	T3N0M0	IIA	Malignant	
E4	44	63	F	Colon	Neuroendocrine carcinoma	T3N0M0	IIA	Malignant	
E5	45	57	M	Small intestine	Carcinoid	T2N0M0	I	Malignant	
E6	46	57	M	Small intestine	Carcinoid	T2N0M0	I	Malignant	
E7	47	70	F	Rectum	Carcinoid	T3N0M0	IIA	Malignant	
E8	48	70	F	Rectum	Carcinoid	T3N0M0	IIA	Malignant	
E9	49	53	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
E10	50	53	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F1	51	57	M	Lung	Atypical carcinoid	T1N0M0	IB	Malignant	
F2	52	57	M	Lung	Atypical carcinoid	T1N0M0	IB	Malignant	
F3	53	60	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F4	54	60	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F5	55	67	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F6	56	67	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F7	57	60	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F8	58	60	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F9	59	53	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
F10	60	53	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G1	61	58	M	Lung	Atypical carcinoid	T2N1M0	IIA	Malignant	
G2	62	58	M	Lung	Atypical carcinoid	T2N1M0	IIA	Malignant	
G3	63	64	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G4	64	64	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G5	65	72	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G6	66	72	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G7	67	75	F	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G8	68	75	F	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G9	69	59	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
G10	70	59	M	Lung	Atypical carcinoid	T2N0M0	IB	Malignant	
H1	71	71	M	Lung	Carcinoid	T2N0M0	IB	Malignant	
H2	72	71	M	Lung	Carcinoid	T2N0M0	IB	Malignant	

H4	74	48	F	Lung	Carcinoid	T2N0M0		IB	Malignant	
H5	75	52	F	Lung	Carcinoid	T2N0M0		IB	Malignant	
H6	76	52	F	Lung	Carcinoid	T2N0M0		IB	Malignant	
H7	77	42	M	Lung	Atypical carcinoid	T2N1M0		IIA	Malignant	
H8	78	42	M	Lung	Atypical carcinoid	T2N1M0		IIA	Malignant	
H9	79	39	M	Lung	Carcinoid	T2N0M0		IB	Malignant	
H10	80	39	M	Lung	Carcinoid	T2N0M0		IB	Malignant	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
 184 Park Drive, Milton Park
 Abingdon OX14 4SE.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

North America 
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

Europe 
 Berenkoog 41,
 1822 BH Alkmaar,
 Netherlands
 T: +31 (0) 72 8080244
 F: +31 (0) 72 8080142

Switzerland 
 Via Lisano 3,
 (CP.683)
 CH-6900
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Germany 
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880