



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

GL631 : Medulloblastoma of cerebellum with cerebellum tissue array, 23 cases/63 cores, replacing BC17012c.

Microarray Panel	Medulloblastoma of cerebellum with cerebellum tissue microarray, containing 20 cases of medulloblastoma, triplicated cores per case, 3 cerebellum tissue, single core per case	
Cores	63	GL631 H&E 
Cases	23	
Row number	7	
Column number	9	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8	9
A	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
B	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
C	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
D	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
E	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
F	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
G	Cer	Cer	Cer	Cer	Cer	Cer	Ceb	Ceb	Ceb

Legend: Ceb - Cerebellum tissue, Cer - Cerebellum
● - Malignant tumor, ● - Normal tissue

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	27	F	Cerebellum	Medulloblastoma				Malignant	
A2	2	27	F	Cerebellum	Medulloblastoma				Malignant	
A3	3	27	F	Cerebellum	Medulloblastoma				Malignant	
A4	4	41	M	Cerebellum	Medulloblastoma				Malignant	
A5	5	41	M	Cerebellum	Medulloblastoma				Malignant	
A6	6	41	M	Cerebellum	Medulloblastoma				Malignant	
A7	7	34	M	Cerebellum	Medulloblastoma				Malignant	
A8	8	34	M	Cerebellum	Medulloblastoma				Malignant	
A9	9	34	M	Cerebellum	Medulloblastoma				Malignant	
B1	10	7	M	Cerebellum	Medulloblastoma				Malignant	
B2	11	7	M	Cerebellum	Medulloblastoma				Malignant	
B3	12	7	M	Cerebellum	Medulloblastoma				Malignant	
B4	13	12	F	Cerebellum	Medulloblastoma				Malignant	
B5	14	12	F	Cerebellum	Medulloblastoma				Malignant	
B6	15	12	F	Cerebellum	Medulloblastoma				Malignant	

B7	16	8	F	Cerebellum	Medulloblastoma					Malignant	
B8	17	8	F	Cerebellum	Medulloblastoma					Malignant	
B9	18	8	F	Cerebellum	Medulloblastoma					Malignant	
C1	19	9	M	Cerebellum	Medulloblastoma					Malignant	
C2	20	9	M	Cerebellum	Medulloblastoma					Malignant	
C3	21	9	M	Cerebellum	Medulloblastoma					Malignant	
C4	22	47	F	Cerebellum	Medulloblastoma					Malignant	
C5	23	47	F	Cerebellum	Medulloblastoma					Malignant	
C6	24	47	F	Cerebellum	Medulloblastoma					Malignant	
C7	25	10	M	Cerebellum	Medulloblastoma					Malignant	
C8	26	10	M	Cerebellum	Medulloblastoma					Malignant	
C9	27	10	M	Cerebellum	Medulloblastoma					Malignant	
D1	28	32	F	Cerebellum	Medulloblastoma					Malignant	
D2	29	32	F	Cerebellum	Medulloblastoma					Malignant	
D3	30	32	F	Cerebellum	Medulloblastoma					Malignant	
D4	31	42	M	Cerebellum	Medulloblastoma					Malignant	
D5	32	42	M	Cerebellum	Medulloblastoma					Malignant	
D6	33	42	M	Cerebellum	Medulloblastoma					Malignant	
D7	34	6	M	Cerebellum	Medulloblastoma					Malignant	
D8	35	6	M	Cerebellum	Medulloblastoma					Malignant	
D9	36	6	M	Cerebellum	Medulloblastoma					Malignant	
E1	37	4	F	Cerebellum	Medulloblastoma					Malignant	
E2	38	4	F	Cerebellum	Medulloblastoma					Malignant	
E3	39	4	F	Cerebellum	Medulloblastoma					Malignant	
E4	40	6	F	Cerebellum	Medulloblastoma					Malignant	
E5	41	6	F	Cerebellum	Medulloblastoma					Malignant	
E6	42	6	F	Cerebellum	Medulloblastoma					Malignant	
E7	43	6	M	Cerebellum	Medulloblastoma					Malignant	
E8	44	6	M	Cerebellum	Medulloblastoma					Malignant	
E9	45	6	M	Cerebellum	Medulloblastoma					Malignant	
F1	46	41	M	Cerebellum	Medulloblastoma					Malignant	
F2	47	41	M	Cerebellum	Medulloblastoma					Malignant	
F3	48	41	M	Cerebellum	Medulloblastoma					Malignant	
F4	49	49	M	Cerebellum	Medulloblastoma					Malignant	
F5	50	49	M	Cerebellum	Medulloblastoma					Malignant	
F6	51	49	M	Cerebellum	Medulloblastoma					Malignant	
F7	52	18	F	Cerebellum	Medulloblastoma					Malignant	
F8	53	18	F	Cerebellum	Medulloblastoma					Malignant	
F9	54	18	F	Cerebellum	Medulloblastoma					Malignant	
G1	55	25	M	Cerebellum	Medulloblastoma					Malignant	
G2	56	25	M	Cerebellum	Medulloblastoma					Malignant	
G3	57	25	M	Cerebellum	Medulloblastoma					Malignant	
G4	58	11	M	Cerebellum	Medulloblastoma					Malignant	
G5	59	11	M	Cerebellum	Medulloblastoma					Malignant	
G6	60	11	M	Cerebellum	Medulloblastoma					Malignant	
G7	61	39	F	Cerebellum	Cerebellum tissue					Normal	
G8	62	47	M	Cerebellum	Cerebellum tissue					Normal	
G9	63	54	F	Cerebellum	Cerebellum tissue					Normal	
-	0										

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

North America 
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

Europe 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

Switzerland 
Via Lisano 3,
(CP.683)
CH-6900
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

Germany 
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880