



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Cat #: HD-201-HR

Description: Human Brain cDNA-Random Primer

30 reactions

Product Description

Human Brain full length tissue cDNA is made of RNA extracted from freshly harvested tissues of single healthy normal adult donor using classical guanidine isothiocyanate phenol: chloroform extraction method. RNA is treated with RNase-free DNase-1 to remove residual DNA. The cDNA is primed with Random Hexamer primer.

Quality Control:

The integrity of each RNA sample as indicated by intact ribosomal RNA is verified by denatured agarose gel electrophoresis. The purity of RNA is assessed by NanoDrop (A260/A280: 1.9-2.1). Residual DNA contamination is tested by PCR.

The synthesized cDNA is also tested as template for PCR amplification of β -actin gene. PCR product of β -actin was visualized on 1% agarose gel.

Applications:

cDNA is ideal for gene expression analysis by PCR, characterization of alternative splicing of mRNA, verification of genetic mutation, gene cloning and target sequencing.

Packing/shipping:

Each cDNA sample is routinely shipped on dry ice in 1.5 ml vials (30ul each) and is enough for 30 reactions (30 PCR amplifications).

MSDS and Certificate of Analysis in PDF files: Contact info@amsbio.com.

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

North America 
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

Europe 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 728080244
F: +31 (0) 728080142

Switzerland 
Via Lisano 3,
(CP.683)
CH-6900
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

Germany 
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880