



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Human Colon cDNA-Oligo-dT

Cat. HD-311

Human Colon full length tissue cDNA is made of RNA extracted from freshly harvested tissues of single healthy normal donor using classical guanidine isothiocyanate phenol:chloroform extraction method. RNA is treated with RNase-free DNase-1 to remove residual DNA. The cDNA is primed with oligo dT primer.

Quality Control: The integrity of each RNA sample as indicated by intact ribosomal RNA is verified by denatured agarose gel electrophoresis. The purity of RNA is assessed by NanoDrop spectrophotometer (A260/A280: 1.9-2.1). Residual DNA contamination is tested by PCR.

The synthesized cDNA is also used as template for PCR amplification of β -actin gene. PCR product of β -actin was visualized on 1% agarose gel.

Applications: cDNA is ideal for gene expression analysis by PCR, characterization of alternative splicing of mRNA, verification of genetic mutation, gene cloning and target sequencing.

Packing/shipping: each cDNA sample is routinely shipped on dry ice in 1.5 ml vials (30ul each) and is enough for 30 reactions (30 PCR amplifications).

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

North America 
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

Europe 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 728080244
F: +31 (0) 728080142

Switzerland 
Via Lisano 3,
(CP.683)
CH-6900
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

Germany 
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880