



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Human Universal Reference cDNA-Oligo dT

Cat# HD-UR-80

Human Universal Reference cDNA-Oligo dT is a full length first strand cDNA synthesized from AMSBIO universal reference total RNA which is isolated from wide variety of male and female major tissues including brain, heart, kidney, liver, lung, skeletal muscles, pancreas, spleen, stomach, intestine, colon, testis, ovary and uterus.

The cDNA is primed with Oligo dT primer and can be used in gene cloning, target sequencing and as a control for gene expression analysis by regular PCR and real time PCR. Each cDNA product is enough for 40 PCR reactions.

Features:

- **Synthesized from high quality intact universal reference total RNA**
- **Consistent control for standard data set comparisons**
- **Prepared from several major tissues for broad gene coverage**
- **Made at industrial scales to minimize lot-to-lot variation.**

Quality Control: The integrity of each RNA sample used for cDNA synthesis is verified by denatured agarose gel electrophoresis. The purity of RNA is assessed by spectrophotometer (A260/A280: 1.9-2.1). Residual DNA contamination is tested by PCR.

The synthesized cDNA is also used as template PCR amplification of β -actin gene. PCR product of β -actin was visualized on 1% agarose gel.

Applications: cDNA is ideal as a control for gene expression analysis by real time and regular PCR, characterization of alternative splicing of mRNA, verification of genetic mutation, gene cloning and target sequencing.

Packing/shipping: each universal reference cDNA sample is routinely shipped on dry ice in 1.5 ml vials and is enough for 40 reactions (40 PCR amplifications).

MSDS and Certificate of Analysis in PDF files: Contact AMSBIO Technical Support at info@amsbio.com

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85