



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Human Tissue Total Protein Panel, any 10 Tissues

Human tissue total protein (protein lysate) panel is your choice of any 10 Tissues (0.1mg each). Tissues are homogenized in protein lysis buffer supplemented with a cocktail of 7 mammalian protease inhibitors to minimize proteolysis, precisely quantified, and stored at -80oC.

Panels:

HT-010 Human Tissue Total Protein Panel 10 x 0.1mg

Total Proteins Available:

HT-201	Human Brain Total Protein
HT-311	Human Colon Total Protein
HT-801	Human Heart Total Protein
HT-901	Human Kidney Total Protein
HT-314	Human Liver Total Protein
HT-601	Human Lung Total Protein
HT-406	Human Ovary Total Protein
HT-313	Human Pancreas Total Protein
HT-413	Human Placenta Total Protein
HT-408	Human Prostate Total Protein
HT-102	Human Skeletal Muscles Total Protein
HT-101	Human Skin Total Protein
HT-306	Human Small Intestine Total Protein
HT-701	Human Spleen Total Protein
HT-302	Human Stomach Total Protein
HT-401	Human Testis Total Protein
HT-312	Human Rectum Total Protein
HT-101	Human Skin Total Protein

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85

Quality Control: Quality of protein as indicated by the absence of smear (no degradation) and sharpness and resolution of the protein bands is verified by denatured SDS-PAGE using Coomassie blue staining. The integrity of protein is tested by immunoblotting using specific beta-actin antibody.

Applications: Total protein is ready for immediate use in Western blotting, immunoprecipitation, SDS-PAGE, isoelectric focusing gels and SDS-capillary electrophoresis, enzymatic activity analysis, gel shifting assay, protein-protein interaction, and tissue specific expression.

Packing/shipping: The extracted protein is provided in vials at a concentration of 5mg/ml and shipped on dry ice. Protein can also be provided in denatured SDS Sample buffer for only Western blotting use

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85