



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

LY481a : Lymph node tissue array, 24 cases/48 cores, replacing LY481

Microarray Panel	Lymph node tissue microarray, containing 21 cases of lymph node tissue, 3 adjacent normal lymph node tissue, duplicated cores per case		
Cores	48	LY481a H&E 	
Cases	24		
Row number	6		
Column number	8		
Core Diameter (mm)	1.5		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.		

	1	2	3	4	5	6	7	8
A	Lym	Lym	Lym	Lym	Lym	Lym	Lym	Lym
B	Lym	Lym	Lym	Lym	Lym	Lym	Lym	Lym
C	Lym	Lym	Lym	Lym	Lym	Lym	Lym	Lym
D	Lym	Lym	Lym	Lym	Lym	Lym	Lym	Lym
E	Lym	Lym	Lym	Lym	Lym	Lym	Lym	Lym
F	Lym	Lym	Lym	Lym	Lym	Lym	Lym	Lym

Legend:

Lym - Lymph node

● - Malignant tumor, ● - NAT, ● - Normal tissue

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	30	M	Lymph node	Lymph node tissue				Normal	
A2	2	30	M	Lymph node	Lymph node tissue				Normal	
A3	3	34	M	Lymph node	Lymph node tissue				Normal	
A4	4	34	M	Lymph node	Lymph node tissue				Normal	
A5	5	30	M	Lymph node	Lymph node tissue				Normal	
A6	6	30	M	Lymph node	Lymph node tissue				Normal	
A7	7	40	M	Lymph node	Lymph node tissue				Normal	
A8	8	40	M	Lymph node	Lymph node tissue				Normal	
B1	9	32	M	Lymph node	Lymph node tissue				Normal	
B2	10	32	M	Lymph node	Lymph node tissue				Normal	
B3	11	35	M	Lymph node	Lymph node tissue				Normal	
B4	12	35	M	Lymph node	Lymph node tissue				Normal	
B5	13	28	M	Lymph node	Lymph node tissue				Normal	
B6	14	28	M	Lymph node	Lymph node tissue				Normal	
B7	15	38	M	Lymph node	Lymph node tissue				Normal	
B8	16	38	M	Lymph node	Lymph node tissue				Normal	

C1	17	30	M	Lymph node	Lymph node tissue				Normal	
C2	18	30	M	Lymph node	Lymph node tissue				Normal	
C3	19	35	M	Lymph node	Lymph node tissue				Normal	
C4	20	35	M	Lymph node	Lymph node tissue				Normal	
C5	21	36	M	Lymph node	Lymph node tissue				Normal	
C6	22	36	M	Lymph node	Lymph node tissue				Normal	
C7	23	31	M	Lymph node	Lymph node tissue				Normal	
C8	24	31	M	Lymph node	Lymph node tissue				Normal	
D1	25	28	M	Lymph node	Lymph node tissue				Normal	
D2	26	28	M	Lymph node	Lymph node tissue				Normal	
D3	27	31	M	Lymph node	Lymph node tissue				Normal	
D4	28	31	M	Lymph node	Lymph node tissue				Normal	
D5	29	33	M	Lymph node	Lymph node tissue				Normal	
D6	30	33	M	Lymph node	Lymph node tissue				Normal	
D7	31	40	M	Lymph node	Lymph node tissue				Normal	
D8	32	40	M	Lymph node	Lymph node tissue				Normal	
E1	33	28	M	Lymph node	Lymph node tissue				Normal	
E2	34	28	M	Lymph node	Lymph node tissue				Normal	
E3	35	32	M	Lymph node	Lymph node tissue				Normal	
E4	36	32	M	Lymph node	Lymph node tissue				Normal	
E5	37	32	M	Lymph node	Lymph node tissue				Normal	
E6	38	32	M	Lymph node	Lymph node tissue				Normal	
E7	39	45	M	Lymph node	Lymph node tissue				Normal	
E8	40	45	M	Lymph node	Lymph node tissue				Normal	
F1	41	35	M	Lymph node	Lymph node tissue				Normal	
F2	42	35	M	Lymph node	Lymph node tissue				Normal	
F3	43	79	F	Lymph node	Adjacent normal lymph node tissue of armpit				NAT	
F4	44	79	F	Lymph node	Adjacent normal lymph node tissue of armpit				NAT	
F5	45	42	M	Lymph node	Adjacent normal lymph node tissue of palate				NAT	
F6	46	42	M	Lymph node	Adjacent normal lymph node tissue of palate				NAT	
F7	47	56	M	Lymph node	Adjacent normal lymph node tissue of mesocolon				NAT	
F8	48	56	M	Lymph node	Adjacent normal lymph node tissue of mesocolon				NAT	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada

1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU

Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World

184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland

Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland

T: +49 (0) 69 779099
F: +49 (0) 69 13376880