



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Mouse Universal Reference cDNA-Random Hexamer

Cat. MD-UR-80-RH

Mouse Universal Reference RNA is composed of highly pure intact total RNA isolated from major male and female C57BL/6 mouse tissues including brain, heart, kidney, liver, lung, skeletal muscles, pancreas, spleen, stomach, intestine, colon, testis, ovary and uterus. Total RNA is extracted using modified guanidine thiocyanate technique which allows the rapid isolation of total RNA including microRNAs. Tissues are cut in very small pieces, snap frozen, grinding to powder, and homogenized in RNA denatured solution. The extracted RNA is treated with DNase, packed at a concentration of 1mg/ml, and stored at -80°C

Features:

- High quality intact RNA without genomic DNA contamination
- Consistent control for standard data set comparisons
- Extracted from large variety of male and female human major tissues for broad gene coverage
- Made at industrial scales to minimize lot-to-lot variation

Quality Control: The integrity of RNA sample as indicated by intact ribosomal RNA is verified by denatured agarose gel electrophoresis. The purity of RNA is assessed by spectrophotometer (A260/A280: 1.9 2.1). Residual DNA contamination is tested by PCR.

Applications: Reference RNA is ideal for microarray analysis, Northern blotting, ribonuclease protection assay, SI nuclease assay, RT PCR/Q PCR/RACE analysis, cDNA synthesis, RNA differential display, microRNA studies, and purification of mRNA for library construction.

Packing/shipping: Total RNA sample is routinely provided in RNase-free water at a concentration of 1 mg/ml and shipped on dry ice.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85