



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Ready-to-Use Mouse BALBC Long Bone Frozen Sections

Catalog #: MF-107-BLC

Quantity: 10 slides

Storage Conditions: Keep at -20°C for short-term storage and at -70°C for long-term storage in a tightly closed box or sealed bag. Sections are good for up to 1 year of the date of purchase if stored properly.

Applications: Mouse frozen tissue sections are ideal for:

- Cellular localization of the protein of interest in mouse tissues by immunohistochemistry.
- Cellular localization of the mRNA of interest in mouse tissues by *in situ* hybridization or *in situ* PCR using both radioactive- and non-radioactive-labeled probes.

Description: Bones are freshly harvested from adult mice, fixed for overnight in 10% Neutral buffered formalin, decalcified to soften bone tissue, cryoprotected in serial sucrose solutions and snapped frozen in O.C.T. freezing medium. Frozen tissues are cut at a thickness of 10 microns and mounted on positively charged slides. Each slide contains 2 sections. All tissue preparation steps from harvesting to embedding in OCT are done in RNase-, DNase-, and protease-free conditions.

Shipping: Slides are shipped for O/N delivery on dry ice in a slide mailer or slide box.

Prior to use: Remove slide box/slide mailer from the freezer (don't open) and warm-up to the room temperature for at least 30 minutes. This step is necessary for prevention of water condensation on the slides. After opening the slide box or slide mailer, it is also recommended to leave slides at room temperature for an additional 10-15 minutes. Restore unused slides in the freezer as described above.

It is recommended to fix sections on slides for at least 10min before use.

Intended use of the products: All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval.

Limited product warranty: AMSBIO has a total quality assurance policy. If any product does not meet the specification stated in the product datasheet will be replaced at no charge. This warranty limits AMSBIO's liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not be liable for any incidental, indirect, special or consequential damages, even if it is aware of the possibility of such damages.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85