



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



MG801b : Meningioma-CNS tumor tissue array, 80 cases/80 cores, replacing MG801a

Microarray Panel	Meningioma-CNS Tumor tissue array, Including multiple histological types and WHO classification (2016 ed) of the meningioma, single core per case	
Cores	80	MG801b H&E 
Cases	80	
Row number	8	
Column number	10	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “*” means invalid core; “-” means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8	9	10
A	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
B	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
C	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
D	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
E	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
F	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
G	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer
H	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Cer	Adr

Legend: Cer - Cerebrum convex/Right parietal lobe

● - Benign tumor, ● - Malignant tumor

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	29	F	Cerebrum convex/Frontal parietal lobe	Meningioma, meningothelial meningioma		1		Benign	
A2	2	39	M	Cerebrum convex	Meningioma, meningothelial meningioma		1		Benign	
A3	3	47	F	Cerebrum convex/Sphenoid bone	Meningioma, meningothelial meningioma		1		Benign	
A4	4	56	M	Cerebrum convex	Meningioma, meningothelial meningioma		1		Benign	
A5	5	67	F	Cerebrum convex/Saddle side	Meningioma, meningothelial meningioma		1		Benign	
A6	6	44	F	Cerebrum convex/Left temporal lobe	Meningioma, meningothelial meningioma		1		Benign	
A7	7	47	F	Cerebrum convex	Meningioma, meningothelial meningioma		1		Benign	
A8	8	67	M	Cerebrum convex/Right frontal	Meningioma, meningothelial meningioma		1		Benign	

A9	9	43	M	Cerebrum convex/Parietal lobe	Meningioma, meningothelial meningioma	1		Benign	
A10	10	45	M	Cerebrum convex/Right parietal lobe	Meningioma, meningothelial meningioma	1		Benign	
B1	11	27	M	Cerebrum convex/Saddle area	Meningioma, meningothelial meningioma	1		Benign	
B2	12	50	F	Cerebrum convex/Left temporal lobe	Meningioma, meningothelial meningioma	1		Benign	
B3	13	52	M	Cerebrum convex/Anterior fossa	Meningioma, meningothelial meningioma	1		Benign	
B4	14	68	F	Cerebrum convex/Left frontal crest	Meningioma, meningothelial meningioma	1		Benign	
B5	15	35	F	Cerebrum convex/Olfactory groove	Meningioma, meningothelial meningioma	1		Benign	
B6	16	56	M	Cerebrum convex	Meningioma, meningothelial meningioma	1		Benign	
B7	17	16	M	Cerebrum convex/Right temporal lobe	Meningioma, meningothelial meningioma	1		Benign	
B8	18	62	M	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
B9	19	35	F	Cerebrum convex/Right frontal	Meningioma, fibrous meningioma	1		Benign	
B10	20	61	F	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
C1	21	52	F	Cerebrum convex/Right cerebellar pontine Angle	Meningioma, fibrous meningioma	1		Benign	
C2	22	27	M	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
C3	23	50	M	Cerebrum convex/Right temporal lobe	Meningioma, fibrous meningioma	1		Benign	
C4	24	37	F	Cerebrum convex/Left occipital lobe	Meningioma, fibrous meningioma	1		Benign	
C5	25	48	F	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
C6	26	40	F	Cerebrum convex/Frontal parietal lobe	Meningioma, fibrous meningioma	1		Benign	
C7	27	38	M	Cerebrum convex/Left parietal lobe	Meningioma, fibrous meningioma	1		Benign	
C8	28	60	F	Cerebrum convex/Right temporal lobe	Meningioma, fibrous meningioma	1		Benign	
C9	29	74	F	Cerebrum convex/Left falx	Meningioma, fibrous meningioma	1		Benign	
C10	30	61	F	Cerebrum convex/Sphenoid ridge	Meningioma, fibrous meningioma	1		Benign	
D1	31	40	F	Cerebrum convex/Right frontal	Meningioma, fibrous meningioma	1		Benign	
D2	32	51	F	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
D3	33	54	M	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
D4	34	62	F	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
D5	35	40	M	Cerebrum convex/Left frontal	Meningioma, fibrous meningioma	1		Benign	
D6	36	58	F	Cerebrum convex/Left temporal lobe	Meningioma, fibrous meningioma	1		Benign	
D7	37	43	M	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
D8	38	67	F	Cerebrum convex	Meningioma, fibrous meningioma	1		Benign	
D9	39	12	M	Cerebrum convex/Left temporal lobe	Meningioma, fibrous meningioma	1		Benign	

D10	40	61	F	Cerebrum convex	Meningioma, fibrous meningioma	1	Benign	
E1	41	75	F	Cerebrum convex/Left parietal lobe	Meningioma, fibrous meningioma	1	Benign	
E2	42	52	F	Cerebrum convex/Left parietal sinus	Meningioma, fibrous meningioma	1	Benign	
E3	43	38	F	Cerebrum convex/Right parietal lobe	Meningioma, fibrous meningioma	1	Benign	
E4	44	48	F	Cerebrum convex/Right parietal lobe	Meningioma, fibrous meningioma	1	Benign	
E5	45	32	F	Cerebrum convex	Meningioma, fibrous meningioma	1	Benign	
E6	46	14	F	Cerebrum convex/Left frontal	Meningioma, fibrous meningioma	1	Benign	
E7	47	60	M	Cerebrum convex	Meningioma, transitional meningioma	1	Benign	
E8	48	66	M	Cerebrum convex	Meningioma, fibrous meningioma	1	Benign	
E9	49	78	M	Cerebrum convex	Meningioma, fibrous meningioma	1	Benign	
E10	50	57	F	Cerebrum convex	Meningioma, fibrous meningioma	1	Benign	
F1	51	39	M	Cerebrum convex/Saddle area	Meningioma, transitional meningioma	1	Benign	
F2	52	59	F	Cerebrum convex/Right sphenoid ridge	Meningioma, fibrous meningioma	1	Benign	
F3	53	40	M	Cerebrum convex/Anterior fossa	Meningioma, fibrous meningioma	1	Benign	
F4	54	72	F	Cerebrum convex	Meningioma, transitional meningioma	1	Benign	
F5	55	40	F	Cerebrum convex/Frontal	Meningioma, fibrous meningioma	1	Benign	
F6	56	36	F	Cerebrum convex/Right ventricle	Meningioma, fibrous meningioma	1	Benign	
F7	57	43	M	Cerebrum convex/Occipital lobe	Meningioma, fibrous meningioma	1	Benign	
F8	58	48	F	Cerebrum convex/Frontal	Meningioma, fibrous meningioma	1	Benign	
F9	59	53	F	Cerebrum convex/Olfactory groove	Meningioma, psammomatous meningioma	1	Benign	
F10	60	66	M	Cerebrum convex	Meningioma, transitional meningioma	1	Benign	
G1	61	51	F	Cerebrum convex	Meningioma, angiomatous meningioma	1	Benign	
G2	62	49	F	Cerebrum convex	Meningioma, microcystic meningioma	1	Benign	
G3	63	64	F	Cerebrum convex/Left temporal lobe	Meningioma, chordoid meningioma	2	Benign	
G4	64	38	M	Cerebrum convex/Anterior fossa	Meningioma, transitional meningioma	1	Benign	
G5	65	74	M	Cerebrum convex/Anterior fossa	Meningioma, transitional meningioma	1	Benign	
G6	66	44	M	Cerebrum convex/Right frontal	Meningioma, transitional meningioma	1	Benign	
G7	67	46	M	Cerebrum convex/Left ventricle	Meningioma, transitional meningioma	1	Benign	
G8	68	57	M	Cerebrum convex/Cerebellar pedicle	Meningioma, atypical meningioma	2	Benign	
G9	69	63	F	Cerebrum convex/Left falx/parafalcine	Meningioma, angiomatous meningioma	1	Benign	
G10	70	68	F	Cerebrum convex/Left frontal	Meningioma, transitional meningioma	1	Benign	
H1	71	65	F	Cerebrum convex/Right parietal lobe	Meningioma, fibrous meningioma	1	Benign	
H2	72	64	F	Cerebrum convex/Right temporal lobe	Meningioma, meningothelial meningioma	1	Benign	

				obe					
H3	73	38	F	Cerebrum convex/Saddle area	Meningioma, papillary meningioma		3		Malignant
H4	74	53	M	Cerebrum convex	Meningioma, transitional meningioma		1		Benign
H5	75	50	F	Cerebrum convex	Meningioma, transitional meningioma		1		Benign
H6	76	62	F	Cerebrum convex/Right frontal	Meningioma, transitional meningioma		1		Benign
H7	77	36	F	Cerebrum convex	Meningioma, transitional meningioma		1		Benign
H8	78	72	M	Cerebrum convex	Meningioma, atypical meningioma		2		Benign
H9	79	57	M	Cerebrum convex/Right temporal lobe	Meningioma, atypical meningioma		2		Benign
H10	80	47	M	Cerebrum convex/Right parietal lobe	Meningioma, atypical meningioma		2		Benign
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC

USA & Canada



1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV

EU



Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd

UK & Rest of the World



184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd

Switzerland



Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV

Deutschland



T: +49 (0) 69 779099
F: +49 (0) 69 13376880

©AMSBIO / AMS Biotechnology (Europe) Ltd. AMSBIO is the global source of all the products listed here.

All trademarks mentioned here are the property of their respective owners.