



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Mouse C57 Brain Region 9 Paraffin Sections

Catalog #: MP-201-09-C57

Quantity : 10 slides

Storage Conditions: Store slides at room temperature in a dry place or at +4C in a sealed bag containing desiccant. Sections are good for 1 year of the date of purchase, if stored properly.

Applications: Mouse paraffin tissue sections ideal for:

- Cellular localization of the protein of interest in mouse embryonic tissues by immunohistochemistry
- Cellular localization of the mRNA of interest in mouse embryonic tissues by *in situ* hybridization or *in situ* PCR using both radioactive- and non-radioactive labelled probes.

Product Description

Mouse tissues are freshly harvested from C57 mice, fixed in 10% neutral buffered formalin and processed for embedding in low temperature melting paraffin. Tissue blocks cut to a thickness of 5microns and mounted on positively charged slides. Each slide contains 2-3 sections depending on the size of sections. All tissue preparation steps from harvesting to embedding in paraffin are done in RNase-, DNase-, and protease-free conditions.

Shipping: Slides are shipped at room temperature in a slide mailer or slide box.

Prior to use: Remove paraffin by dipping slides in xylene or any substitutive agent (3 changes, 5 minutes each), rehydrate in serial ethanol dilutions, and finally wash with PBS. For *in situ* hybridization applications, make alcohol solutions and PBS in RNase-free water.

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880