



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

*Data Sheet***Ready-to-use Mouse-Balbc Brain Sagittal Paraffin Sections****Catalog #:** MP-201-SS-BLC**Quantity:** 10 slides**Storage Conditions:** Store slides at room temperature in a dry place or at 4°C in a sealed bag containing desiccant. Sections are good for 1 year of the date of purchase if stored properly.**Applications:** Paraffin tissue sections are ideal for:

- Cellular localization of protein of interest in mouse tissues by immunohistochemistry.
- Cellular localization of mRNA of interest in mouse tissues by *in situ* hybridization or *in situ* PCR using both radioactive- and non-radioactive-labeled probes.

Description: Brain is freshly harvested from adult Balbc mouse, fixed in 10% buffered neutral formalin and processed for embedding in low temperature melting paraffin. Tissue blocks are cut at a thickness of 5 microns and mounted on positively charge slides. All tissue preparation steps from harvesting to embedding in paraffin are done in RNase-, DNase-, and protease-free conditions.**Shipping:** Slides are shipped at room temperature in a slide mailer or slide box.**Prior to use:** Remove paraffin by dipping slides in xylene or any substitutive agent (3 changes 5 minute each), rehydrate in serial ethanol dilutions, and finally wash with PBS. For *in situ* hybridization use, make alcohol solutions and PBS in RNase-free water.**Intended use of the products:** All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. Our products may not resell, modified for resale or used to manufacture commercial products without written approval.**Limited product warranty:** AMSBIO has a total quality assurance policy. If any product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not be liable for any incidental, indirect, special or consequential damages, even if it is aware of the possibility of such damages.AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482



North America
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218



Germany
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880



Switzerland
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85