



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse CD1 Embryo Total RNA Panel - Set of 5 Embryos

Catalog #: MR-104-005

Quantity: 5x0.025mg

Storage Conditions: Store at -80°C. Expiry date is one year of the date of purchase if stored correctly.

Applications: RNA is ideal for Northern blotting, ribonuclease protection assay, SI nuclease assay, RT-PCR analysis, rapid amplification of cDNA ends (RACE) and purification of mRNA for library construction.

Quality Control: The integrity of each RNA sample as indicated by intact 28S and 18S ribosomal RNA is verified by denaturing agarose gel electrophoresis. The purity of RNA is assessed by spectrophotometer ($A_{260}/A_{280} > 1.9$).

Description: Total RNA is isolated from 5 mouse tissues (Balb/C strain) selected by the customer using a modified guanidine thiocyanate technique. RNA samples are treated with RNase-free DNase to remove residual DNA, suspended in water, precisely quantified, and stored at -80°C.

List of Stages Available:

Mouse CD1 Embryo-E10 Total RNA
Mouse CD1 Embryo-E11 Total RNA
Mouse CD1 Embryo-E12 Total RNA
Mouse CD1 Embryo-E13 Total RNA
Mouse CD1 Embryo-E14 Total RNA
Mouse CD1 Embryo-E15 Total RNA
Mouse CD1 Embryo-E16 Total RNA
Mouse CD1 Embryo-E17 Total RNA
Mouse CD1 Embryo-E18 Total RNA

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85