



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Rat WS Lung Frozen Sections

Catalog #: RF-601-WS

Quantity: 10 slides

Storage Conditions: Keep at -20oC for short-term storage and at -70oC for long-term storage in a tightly closed box or sealed bag. Sections are good for up to 1 year of the date of purchase if stored properly.

Applications: Rat frozen tissue sections are ideal for:

- Cellular localization of the protein of interest in tissues by immunohistochemistry.
- Cellular localization of the mRNA of interest in tissues by *in situ* hybridization or *in situ* PCR using both radioactive- and non-radioactive-labeled probes.

Description: Tissues are freshly harvested from Wistar rat and snapped frozen in O.C.T. freezing medium. Frozen tissues are cut at a thickness of 7-10 microns and mounted on slides treated for adherence. Each slide contains 1 to 4 sections depending on the size of the frozen tissue. All tissue preparation steps from harvesting to embedding in OCT medium are done in RNase-, DNase-, and protease-free conditions.

Shipping: Slides are shipped for O/N delivery on dry ice in a slide mailer or slide box.

Prior to use: Remove slide box/slide mailer from the freezer (don't open) and warm-up to the room temperature for at least 30 minutes. This step is necessary for prevention of water condensation on the slides. After opening the slide box or slide mailer, it is also recommended to leave slides at room temperature for an additional 10-15 minutes. Restore unused slides in the freezer as described above.

Tissues are not yet fixed and need to be fixed in one of the following fixatives depending on your applications: ice-cold acetone, ethanol, 4%paraformaldehyde, 10% buffered neutral formalin, or any other fixatives recommended for frozen tissue sections.

Intended use of the products: All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval from AMSBIO.
Limited product warranty: AMSBIO has a total quality assurance policy. If any product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not liable for any incidental, indirect, special or consequential damages, even if it is aware of the possibility of such damages.

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
 184 Park Drive, Milton Park
 Abingdon OX14 4SE.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

North America 
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

Europe 
 Berenkoog 41,
 1822 BH Alkmaar,
 Netherlands
 T: +31 (0)728080244
 F: +31 (0)728080142

Switzerland 
 Via Lisano 3,
 (CP.683)
 CH-6900
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Germany 
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880