



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

SS1001 : Synovial sarcoma tissue array, including TNM, clinical stage and pathology grade, 50 cases/100 cores

Microarray Panel	Synovial sarcoma tissue microarray, duplicate core per case		
Cores	100		
Cases	50		
Row number	10		
Column number	10		
Core Diameter (mm)	1		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "*" means invalid core; "-" means no applicable or negative in IHC markers.		

	1	2	3	4	5	6	7	8	9	10
A	Art	Art	Art	Art	Art	Art	Fib	Fib	Fib	Fib
B	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten
C	Ten	Ten	Bon	Bon	Ner	Ner	Blo	Blo	Fib	Fib
D	Fib	Fib	Fib	Fib	Fib	Fib	Adi	Adi	Str	Str
E	Smo	Smo	Str	Str	Ten	Ten	Ten	Ten	Ten	Ten
F	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten
G	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten
H	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten
I	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten	Ten
J	Ten	Ten	Sof	Sof	Sof	Sof	Sof	Sof	Sof	Adr

Legend:Adi - Adipose tissue, Art - Articulation, Blo - Blood vessel, Bon - Bone, Fib - Fibrous tissue, Ner - Nerve, Smo - Smooth muscle, Sof - Soft tissues, Str - Striated muscle, Ten - Tendon sheath

● - Malignant tumor, ○ - Malignant tumor (stage IA), ○ - Malignant tumor (stage IB), ○ - Malignant tumor (stage IIA), ○ - Malignant tumor (stage IIB), ○ - Malignant tumor (stage IIIA), ○ - Malignant tumor (stage IIIB), ○ - Malignant tumor (stage IVA)

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	53	F	Articulation	Synovial sarcoma of lefe articular genu	T2N0M0 G3	-	IIIB	malignant	
A2	2	53	F	Articulation	Synovial sarcoma of lefe articular genu	T2N0M0 G3	-	IIIB	malignant	
A3	3	36	F	Articulation	Synovial sarcoma of lefe sacroiliac joint	T1aN0M0 G1	-	IA	malignant	
A4	4	36	F	Articulation	Synovial sarcoma of lefe sacroiliac joint	T1aN0M0 G1	-	IA	malignant	
A5	5	23	F	Articulation	Synovial sarcoma of right articular genu	T2bN0M0 G3	-	IIIB	malignant	
A6	6	23	F	Articulation	Synovial sarcoma of right articular genu	T2bN0M0 G3	-	IIIB	malignant	
A7	7	14	F	Fibrous tissue	Synovial sarcoma of lefe articular genu	T1aN0M0 G3	-	IIIA	malignant	
A8	8	14	F	Fibrous tissue	Synovial sarcoma of lefe articular genu	T1aN0M0 G3	-	IIIA	malignant	
A9	9	64	M	Fibrous tissue	Synovial sarcoma of lefe articular genu	T2N0M0 G2	-	IIB	malignant	
A10	10	64	M	Fibrous tissue	Synovial sarcoma of lefe articular genu	T2N0M0 G2	-	IIB	malignant	
B1	11	37	F	Tendon sheath	Synovial sarcoma of right wrist joint (sparse)	T1aN0M0 G1	-	IA	malignant	
B2	12	37	F	Tendon sheath	Synovial sarcoma of right wrist joint (sparse)	T1aN0M0 G1	-	IA	malignant	
B3	13	23	M	Tendon sheath	Synovial sarcoma of right articular genu	T2aN0M0 G1	-	IB	malignant	
B4	14	23	M	Tendon sheath	Synovial sarcoma of right articular genu	T2aN0M0 G1	-	IB	malignant	
B5	15	46	M	Tendon sheath	Synovial sarcoma of left articular genu	T1aN0M0 G1	-	IA	malignant	
B6	16	46	M	Tendon sheath	Synovial sarcoma of left articular genu	T1aN0M0 G1	-	IA	malignant	
B7	17	58	M	Tendon sheath	Synovial sarcoma of left articular genu	T2N0M0 G3	-	IIIB	malignant	
B8	18	58	M	Tendon sheath	Synovial sarcoma of left articular genu	T2N0M0 G3	-	IIIB	malignant	
B9	19	28	M	Tendon sheath	Synovial sarcoma of right articular genu	T1aN0M0 G2	-	IIA	malignant	
B10	20	28	M	Tendon sheath	Synovial sarcoma of right articular genu	T1aN0M0 G2	-	IIA	malignant	

C1	21	30	F	Tendon sheath	Synovial sarcoma of left articular genu	T1aN0M0 G2	-	IIA	malignant
C2	22	30	F	Tendon sheath	Synovial sarcoma of left articular genu	T1aN0M0 G2	-	IIA	malignant
C3	23	35	M	Bone	Synovial sarcoma of right articular genu	T2N0M0 G1	-	IB	malignant
C4	24	35	M	Bone	Synovial sarcoma of right articular genu	T2N0M0 G1	-	IB	malignant
C5	25	46	M	Nerve	Synovial sarcoma of right iliac region	T2N1M0 G3	-	IVA	malignant
C6	26	46	M	Nerve	Synovial sarcoma of right iliac region	T2N1M0 G3	-	IVA	malignant
C7	27	43	M	Blood vessel	Synovial sarcoma of right thigh	T1bN0M0 G3	-	IIIA	malignant
C8	28	43	M	Blood vessel	Synovial sarcoma of right thigh	T1bN0M0 G3	-	IIIA	malignant
C9	29	48	M	Fibrous tissue	Synovial sarcoma of left thigh	T1bN0M0 G3	-	IIIA	malignant
C10	30	48	M	Fibrous tissue	Synovial sarcoma of left thigh	T1bN0M0 G3	-	IIIA	malignant
D1	31	40	M	Fibrous tissue	Synovial sarcoma of right thigh	T1bN0M0 G3	-	IIIA	malignant
D2	32	40	M	Fibrous tissue	Synovial sarcoma of right thigh	T1bN0M0 G3	-	IIIA	malignant
D3	33	34	M	Fibrous tissue	Synovial sarcoma of left thigh	T2bN0M0 G3	-	IIIB	malignant
D4	34	34	M	Fibrous tissue	Synovial sarcoma of left thigh	T2bN0M0 G3	-	IIIB	malignant
D5	35	54	F	Fibrous tissue	Synovial sarcoma of left leg	T2bN0M0 G3	-	IIIB	malignant
D6	36	54	F	Fibrous tissue	Synovial sarcoma of left leg	T2bN0M0 G3	-	IIIB	malignant
D7	37	32	M	Adipose tissue	Synovial sarcoma of right thigh (skeletal muscle tissue)	T2bN0M0 G2	-	IIB	malignant
D8	38	32	M	Adipose tissue	Synovial sarcoma of right thigh (sparse)	T2bN0M0 G2	-	IIB	malignant
D9	39	47	M	Striated muscle	Synovial sarcoma of left leg	T1bN0M0 G3	-	IIIA	malignant
D10	40	47	M	Striated muscle	Synovial sarcoma of left leg	T1bN0M0 G3	-	IIIA	malignant
E1	41	39	M	Smooth muscle	Synovial sarcoma of right thigh	T2bN0M0 G2	-	IIB	malignant
E2	42	39	M	Smooth muscle	Synovial sarcoma of right thigh	T2bN0M0 G2	-	IIB	malignant
E3	43	16	M	Striated muscle	Synovial sarcoma of right lumbosacral area	T2bN0M0 G2	-	IIB	malignant
E4	44	16	M	Striated muscle	Synovial sarcoma of right lumbosacral area	T2bN0M0 G2	-	IIB	malignant
E5	45	25	M	Tendon sheath	Synovial sarcoma of right knee	T3bN0M0 G3	-	IIIB	malignant
E6	46	25	M	Tendon sheath	Synovial sarcoma of right knee	T3bN0M0 G3	-	IIIB	malignant
E7	47	80	F	Tendon sheath	Synovial sarcoma of left croup	T2aN0M0 G3	-	IIIB	malignant
E8	48	80	F	Tendon sheath	Synovial sarcoma of left croup	T2aN0M0 G3	-	IIIB	malignant
E9	49	32	F	Tendon sheath	Synovial sarcoma of left thigh	T1aN0M0 G3	-	IIIA	malignant
E10	50	32	F	Tendon sheath	Synovial sarcoma of left thigh	T1aN0M0 G3	-	IIIA	malignant
F1	51	24	M	Tendon sheath	Synovial sarcoma of right shoulder	T2N0M0 G3	-	IIIB	malignant
F2	52	24	M	Tendon sheath	Synovial sarcoma of right shoulder	T2N0M0 G3	-	IIIB	malignant
F3	53	33	F	Tendon sheath	Synovial sarcoma of left tight (fibrous tissue and blood vessel)	T1bN0M0 G1	-	IA	malignant
F4	54	33	F	Tendon sheath	Synovial sarcoma of left tight (fibrous tissue and blood vessel)	T1bN0M0 G1	-	IA	malignant
F5	55	50	M	Tendon sheath	Synovial sarcoma of right planta	T1aN0M0 G1	-	IA	malignant
F6	56	50	M	Tendon sheath	Synovial sarcoma of right planta	T1aN0M0 G1	-	IA	malignant
F7	57	51	F	Tendon sheath	Synovial sarcoma with necrosis of left thigh	T2bN0M0 G2	-	IIB	malignant
F8	58	51	F	Tendon sheath	Synovial sarcoma of left thigh (sparse)	T2bN0M0 G2	-	IIB	malignant
F9	59	50	M	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G2	-	IIB	malignant
F10	60	50	M	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G2	-	IIB	malignant
G1	61	49	F	Tendon sheath	Synovial sarcoma of right thigh	T1bN0M0 G2	-	IIA	malignant
G2	62	49	F	Tendon sheath	Synovial sarcoma of right thigh	T1bN0M0 G2	-	IIA	malignant
G3	63	25	F	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G2	-	IIB	malignant
G4	64	25	F	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G2	-	IIB	malignant
G5	65	21	F	Tendon sheath	Synovial sarcoma of right elbow	T2N0M0 G2	-	IIIB	malignant
G6	66	21	F	Tendon sheath	Synovial sarcoma of right elbow	T2N0M0 G2	-	IIIB	malignant
G7	67	27	M	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G2	-	IIB	malignant
G8	68	27	M	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G2	-	IIB	malignant
G9	69	37	M	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G3	-	IIIB	malignant
G10	70	37	M	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G3	-	IIIB	malignant
H1	71	76	M	Tendon sheath	Synovial sarcoma of left ankle mortise	T1aN0M0 G1	-	IA	malignant
H2	72	76	M	Tendon sheath	Synovial sarcoma of left ankle mortise	T1aN0M0 G1	-	IA	malignant
H3	73	40	F	Tendon sheath	Synovial sarcoma of hibateral thigh	T2N0M0 G3	-	IIIB	malignant
H4	74	40	F	Tendon sheath	Synovial sarcoma of hibateral thigh	T2N0M0 G3	-	IIIB	malignant
H5	75	47	M	Tendon sheath	Synovial sarcoma with necrosis of os pubis (sparse)	T1bN0M0 G1	-	IA	malignant
H6	76	47	M	Tendon sheath	Synovial sarcoma with necrosis of os pubis	T1bN0M0 G1	-	IA	malignant
H7	77	44	F	Tendon sheath	Synovial sarcoma of lefe articular genu	T1aN0M0 G1	-	IA	malignant
H8	78	44	F	Tendon sheath	Synovial sarcoma of lefe articular genu	T1aN0M0 G1	-	IA	malignant
H9	79	24	M	Tendon sheath	Synovial sarcoma of right thigh	T2bN0M0 G2	-	IIB	malignant
H10	80	24	M	Tendon sheath	Synovial sarcoma of right thigh	T2bN0M0 G2	-	IIB	malignant
I1	81	63	M	Tendon sheath	Synovial sarcoma of scapular region	T2bN0M0 G3	-	IIIB	malignant
I2	82	63	M	Tendon sheath	Synovial sarcoma of scapular region	T2bN0M0 G3	-	IIIB	malignant
I3	83	40	F	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G1	-	IB	malignant
I4	84	40	F	Tendon sheath	Synovial sarcoma of left thigh	T2bN0M0 G1	-	IB	malignant
I5	85	23	M	Tendon sheath	Synovial sarcoma of right elbow	T2aN0M0 G2	-	IIB	malignant
I6	86	23	M	Tendon sheath	Synovial sarcoma of right elbow	T2aN0M0 G2	-	IIB	malignant

I7	87	56	M	Tendon sheath	Synovial sarcoma of right thigh	T2aN0M0 G3	-	IIIB	malignant	
I8	88	56	M	Tendon sheath	Synovial sarcoma of right thigh	T2aN0M0 G3	-	IIIB	malignant	
I9	89	69	F	Tendon sheath	Synovial sarcoma of left leg	T1bN0M0 G2	-	IIA	malignant	
I10	90	69	F	Tendon sheath	Synovial sarcoma of left leg	T1bN0M0 G2	-	IIA	malignant	
J1	91	49	M	Tendon sheath	Synovial sarcoma of right thigh	T1aN0M0 G1	-	IA	malignant	
J2	92	49	M	Tendon sheath	Synovial sarcoma of right thigh	T1aN0M0 G1	-	IA	malignant	
J3	93	67	F	Soft tissues	Synovial sarcoma of left leg	T2bN0M0 G1	-	IB	malignant	
J4	94	67	F	Soft tissues	Synovial sarcoma of left leg	T2bN0M0 G1	-	IB	malignant	
J5	95	64	F	Soft tissues	Synovial sarcoma of right thigh	T2N0M0 G2	-	IIB	malignant	
J6	96	64	F	Soft tissues	Synovial sarcoma of right thigh	T2N0M0 G2	-	IIB	malignant	
J7	97	31	M	Soft tissues	Synovial sarcoma of right leg	T2bN0M0 G2	-	IIA	malignant	
J8	98	31	M	Soft tissues	Synovial sarcoma of right leg	T2bN0M0 G2	-	IIA	malignant	
J9	99	38	M	Soft tissues	Synovial sarcoma of hip (sparse)	T2bN0M0 G3	-	IIIB	malignant	
J10	100	38	M	Soft tissues	Synovial sarcoma of hip (chronic inflammation with necrosis of fibrous tissue)	T2bN0M0 G3	-	IIIB	malignant	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880