



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0


F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ST801b : Stomach cancer tissue array with matched adjacent normal stomach tissue, including pathology grade, TNM and clinical stage, 40 cases/80 cores, replacing ST801a

Microarray Panel	Gastric adenocarcinoma with matched adjacent normal tissue microarray, containing 36 cases of adenocarcinoma, 2 mucinous adenocarcinoma, 1 each of signet ring cell carcinoma and carcinoid, duplicate cores per case	
Cores	80	
Cases	40	
Row number	8	
Column number	10	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, "" means invalid core; "-" means no applicable or negative in IHC markers.	

	1	2	3	4	5	6	7	8	9	10
A	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
B	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
C	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
D	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
E	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
F	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
G	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto
H	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto	Sto

Legend:Sto - Stomach

● - Malignant tumor,
 ● - Malignant tumor (stage IB),
 ● - Malignant tumor (stage IIA),
 ● - Malignant tumor (stage IIB),
 ● - Malignant tumor (stage IIIA),
 ● - Malignant tumor (stage IIIB),
 ● - NAT

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	HER2	Image
A1	1	78	F	Stomach	Adenocarcinoma	T3N0M0	1	IIA	Malignant	0	
A2	2	78	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
A3	3	67	M	Stomach	Adenocarcinoma	T2N0M0	1	IB	Malignant	0	
A4	4	67	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
A5	5	59	F	Stomach	Adenocarcinoma	T3N0M0	1	IIA	Malignant	0	
A6	6	59	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
A7	7	50	M	Stomach	Mucinous adenocarcinoma	T3N0M0	1	IIA	Malignant	0	
A8	8	50	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
A9	9	51	F	Stomach	Adenocarcinoma	T3N0M0	1	IIA	Malignant	1+	
A10	10	51	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
B1	11	68	F	Stomach	Adenocarcinoma	T2N0M0	2	IB	Malignant	0	
B2	12	68	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
B3	13	48	M	Stomach	Adenocarcinoma	T2N0M0	1	IB	Malignant	2+	
B4	14	48	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
B5	15	56	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
B6	16	56	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
B7	17	57	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
B8	18	57	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
B9	19	74	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
B10	20	74	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
C1	21	52	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
C2	22	52	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
C3	23	68	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
C4	24	68	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	

C5	25	51	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	3+	
C6	26	51	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
C7	27	55	F	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
C8	28	55	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
C9	29	53	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	3+	
C10	30	53	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
D1	31	54	M	Stomach	Adenocarcinoma	T2N0M0	2	IB	Malignant	3+	
D2	32	54	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
D3	33	69	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
D4	34	69	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
D5	35	72	M	Stomach	Adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
D6	36	72	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
D7	37	56	M	Stomach	Adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
D8	38	56	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
D9	39	53	F	Stomach	Adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
D10	40	53	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
E1	41	61	M	Stomach	Adenocarcinoma	T3N1M0	2	IIB	Malignant	0	
E2	42	61	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
E3	43	49	M	Stomach	Adenocarcinoma	T3N0M0	2	IIA	Malignant	0	
E4	44	49	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
E5	45	53	M	Stomach	Adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
E6	46	53	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
E7	47	54	M	Stomach	Adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
E8	48	54	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
E9	49	50	M	Stomach	Adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
E10	50	50	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
F1	51	63	F	Stomach	Adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
F2	52	63	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
F3	53	68	F	Stomach	Adenocarcinoma	T2N0M0	3	IB	Malignant	0	
F4	54	68	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
F5	55	71	M	Stomach	Adenocarcinoma	T3N1M0	3	IIB	Malignant	0	
F6	56	71	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
F7	57	61	M	Stomach	Adenocarcinoma	T2N0M0	3	IB	Malignant	0	
F8	58	61	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
F9	59	47	M	Stomach	Adenocarcinoma	T2N1M0	3	IIA	Malignant	0	
F10	60	47	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
G1	61	58	M	Stomach	Mucinous adenocarcinoma	T3N0M0	3	IIA	Malignant	0	
G2	62	58	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
G3	63	63	F	Stomach	Adenocarcinoma	T2N0M0	3	IB	Malignant	0	
G4	64	63	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
G5	65	66	M	Stomach	Adenocarcinoma	T2N0M0	3	IB	Malignant	0	
G6	66	66	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
G7	67	52	M	Stomach	Adenocarcinoma	T2N0M0	3	IB	Malignant	0	
G8	68	52	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
G9	69	54	F	Stomach	Adenocarcinoma	T4N1M0	3	IIIA	Malignant	0	
G10	70	54	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
H1	71	51	F	Stomach	Adenocarcinoma	T3N1M0	3	IIB	Malignant	0	
H2	72	51	F	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
H3	73	64	M	Stomach	Signet-ring cell carcinoma	T3N0M0	-	IIA	Malignant	0	
H4	74	64	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
H5	75	73	M	Stomach	Adenocarcinoma	T2N0M0	3	IB	Malignant	0	
H6	76	73	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
H7	77	62	M	Stomach	Carcinoid	T3N1M0	-	IIB	Malignant	0	
H8	78	62	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
H9	79	50	M	Stomach	Adenocarcinoma	T4N0M0	-	IIIB	Malignant	0	
H10	80	50	M	Stomach	Adjacent normal stomach tissue	-	-	-	NAT	0	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)	-	-	-	Malignant		

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880