



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet

Product Name: Rectum Tumor Tissue Array - 64 Different Rectum tumors.

Catalog No.: T8235723-5

Lot No.: A703136

Species: Human

Tissue Type: Tumor

Tissue Array Diagram:

	11	10	9	8	7	6	5	4	3	2	1	
F	O	O	O	O	O	O	O	O	O	O	O	
E	O	O	O	O	O	O	O	O	O	O	O	
D	O	O	O	O	O	O	O	O	O	O	O	
C	O	O	O	O	O	O	O	O	O	O	O	
B	O	O	O	O	O	O	O	O	O	O	O	
A	O	O	O	O	O	O	O	O	O	O	O	

Components:

1. 5 Tissue Array Slides per package
2. Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: 

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85

Data Sheet

Product Name: Tissue Array

Catalog No.: T8235723-5

Lot No.: A703136

Donor Information

Donor Information									
Position	Age/Sex		Organ	Pathological Diagnosis	Tumor History	Tumor Size (cm)	Differentiation	TNM	Tissue ID
A1	29	F	Placenta	Human Normal Placenta					HT00031
A2				Paraffin					
A3	54	F	Rectum	Adenocarcinoma	3 months	6x5x3	Moderately	T3N1M1	HT00564
A4	43	M	Rectum	Adenocarcinoma	1 months	4x4x3	Moderately	T2N1M0	HT00565
A5	49	M	Rectum	Adenocarcinoma	20 days	4x3x3	Moderately	T2N1M0	HT00566
A6	49	M	Rectum	Adenocarcinoma	20 days	5x4x3	Moderately-Well	T3N2M1	HT00567
A7	38	F	Rectum	Adenocarcinoma	20 days	4.5x4x4	Moderately-Well	T3N3M0	HT00568
A8	50	F	Rectum	Adenocarcinoma, ulcer	4 months	5x7	Well	T3N3M0	HT00569
A9	64	M	Rectum	Adenocarcinoma, ulcer	3 months	3x3	Well	T2N1M0	HT00570
A10	60	F	Rectum	Adenocarcinoma	1 years	6.5x6x1	N/A	T3N3M1	HT00571
A11	62	F	Rectum	Adenocarcinoma	3 months+	diameter 6cm	Moderately	T3N2M0	HT00572
B1	46	M	Rectum	Adenocarcinoma, cauliflower like	1 years	diameter 2.5cm	Moderately	T2N0M0	HT00573
B2	51	F	Rectum	Adenocarcinoma, ulcer	2 months+	diameter 5cm	Moderately	T3N1M0	HT00574
B3	66	M	Rectum	Adenocarcinoma, ulcer	3 months	6.5x6x3	Moderately	T3N1M0	HT00575
B4	32	M	Rectum	Adenocarcinoma, ulcer	1 years	6x5	Well	T3N2M0	HT00576
B5	58	F	Rectum	Adenocarcinoma, ulcer	6 months	4.5x3	Moderately	T2N2M0	HT00577
B6	58	M	Rectum	Adenocarcinoma	30 days	4.5x4x3	Well	T2N1M0	HT00578
B7	72	F	Rectum	Squamous cell carcinoma	15 days	3x5	Moderately	T3N0M0	HT00579
B8	43	M	Rectum	Adenocarcinoma	15 days	3.5x3x1.2	Moderately-Well	T2N2M0	HT00580
B9	36	M	Rectum	Adenocarcinoma	10 days	5x3.5x2	Well	T3N2M0	HT00581
B10	59	M	Rectum	Adenocarcinoma	20 days	5x4x3	Moderately-Well	T3N4M1	HT00582
B11	48	M	Rectum	Adenocarcinoma	1 months	4x4x2	Moderately-Poorly	T2N0M0	HT00583
C1	Adult	M	Rectum	Adenocarcinoma	1 months	3x3x1	Moderately-Well	T2N0M0	HT00584
C2	42	M	Rectum	Adenocarcinoma	1 months	4x3.5x1	Moderately-Well	T2N0M0	HT00585
C3	40	F	Rectum	Adenocarcinoma	3 months	2.5x2x0.8	Moderately-Poorly	T2N2M0	HT00586
C4	51	M	Rectum	Adenocarcinoma	15 days	3x3x2	Moderately-Well	T2N0M0	HT00587
C5	50	M	Rectum	Adenocarcinoma	40 days	4x3x3	Moderately	T2N0M0	HT00588
C6	50	F	Rectum	Adenocarcinoma	16 days	2.7x2.5x0.7	Moderately	T2N0M0	HT00589
C7	46	M	Rectum	Adenocarcinoma	24 days	2.8x2.5x0.8	Well	T2N0M0	HT00590
C8	58	M	Rectum	Adenocarcinoma	20 days	5x4x3	Moderately-Well	T3N0M0	HT00591
C9	63	M	Rectum	Adenocarcinoma	1 months	4x3x2.5	Moderately	T2N0M0	HT00592
C10	65	M	Rectum	Adenocarcinoma	50 days	2.7x2.2x1	Moderately-Well	T2N0M0	HT00593
C11	56	M	Rectum	Adenocarcinoma	10 days	2.5x1.5x0.7	Well	T2N0M0	HT00594
D1	46	F	Rectum	Adenocarcinoma	1.5 months	5.5x5x3.5	Poorly	T2N0M0	HT00595

D2	47	M	Rectum	Adenocarcinoma	15 days	2x1.5x1	Well	T ₂ N ₀ M ₀	HT00596
D3	55	M	Rectum	Adenocarcinoma	16 days	2.5x2.5x2	Moderately	T ₂ N ₀ M ₀	HT00597
D4	79	F	Rectum	Adenocarcinoma, ulcer	1 months	5.5x3	Well	T ₃ N ₀ M ₁	HT00598
D5	63	M	Rectum	Adenocarcinoma	20 days	5x4x3	Moderately	T ₃ N ₂ M ₀	HT00599
D6	70	F	Rectum	Adenocarcinoma, mucous	5 months+	4.5X3X1	Well	T ₂ N ₀ M ₀	HT00600
D7	57	F	Rectum	Adenocarcinoma, cauliflower like	2 weeks+	2.5x2x1	Moderately	T ₂ N ₀ M ₀	HT00601
D8	83	M	Rectum	Adenocarcinoma	2 years+	5.5x2.5x2	Moderately	T ₃ N ₁ M ₀	HT00602
D9	70	M	Rectum	Adenosquamous	1 years	3x4	Moderately	T ₂ N ₀ M ₀	HT00603
D10	60	M	Rectum	Adenocarcinoma	3 months	Diameter 3.5cm	N/A	T ₂ N ₀ M ₀	HT00604
D11	23	M	Rectum	Adenocarcinoma	6 months	6x4.5x1.5	Well	T ₃ N ₀ M ₀	HT00605
E1	67	M	Rectum	Adenocarcinoma, ulcer	16 days	5x3.2x1.5	Moderately	T ₃ N ₀ M ₀	HT00606
E2	59	M	Rectum	Adenocarcinoma	8 months	5X4X1.2	Moderately-Well	T ₃ N ₀ M ₀	HT00607
E3	57	M	Rectum	Adenocarcinoma, ulcer	10 days	3x1.5x1	Well	T ₂ N ₀ M ₀	HT00608
E4	62	F	Rectum	Adenocarcinoma, ulcer	6 months	Diameter 3cm	Moderately	T ₂ N ₀ M ₀	HT00609
E5	45	M	Rectum	Adenocarcinoma, ulcer	3 weeks	5.5x4x2	Well	T ₃ N ₀ M ₀	HT00610
E6	60	M	Rectum	Adenocarcinoma, ulcer	1 months	3x2.5x1	Moderately	T ₂ N ₀ M ₀	HT00611
E7	57	M	Rectum	Adenocarcinoma, ulcer	10 months	4.8x4x1.5	Moderately	T ₃ N ₀ M ₀	HT00612
E8	45	M	Rectum	Adenocarcinoma, ulcer	1 months	3x4	Well	T ₂ N ₀ M ₀	HT00613
E9	73	M	Rectum	Adenocarcinoma, ulcer	2 months	6x7	Well	T ₃ N ₀ M ₀	HT00614
E10	64	F	Rectum	Adenocarcinoma, ulcer	9 months	8x3	Well	T ₃ N ₀ M ₀	HT00615
E11	71	F	Rectum	Adenocarcinoma, ulcer	7 days	3x3x2	Well	T ₂ N ₀ M ₀	HT00616
F1	49	M	Rectum	Adenocarcinoma, ulcer & mucous	3 months	4x4	Well	T ₂ N ₀ M ₀	HT00617
F2	43	M	Rectum	Adenocarcinoma, ulcer	2 months	2.5x1.5x1	Poorly	T ₂ N ₀ M ₀	HT00618
F3	57	F	Rectum	Adenocarcinoma, ulcer	2 months	5x3.5x1	Well	T ₃ N ₀ M ₀	HT00619
F4	54	M	Rectum	Adenocarcinoma, ulcer	7 months	1.5x1x1	Moderately	T ₂ N ₀ M ₀	HT00620
F5	63	M	Rectum	Adenocarcinoma, ulcer	2 years	3x3x2	Moderately	T ₂ N ₀ M ₀	HT00621
F6	64	F	Rectum	Adenocarcinoma, ulcer	2 months	Diameter 4cm	Moderately	T ₂ N ₀ M ₀	HT00622
F7	74	M	Rectum	Adenocarcinoma, ulcer	6 months	1.5x1x0.8	Poorly	T ₂ N ₀ M ₀	HT00623
F8	61	M	Rectum	Adenocarcinoma, ulcer	3 months	5x4x2	Moderately	T ₃ N ₀ M ₀	HT00624
F9	40	F	Rectum	Adenocarcinoma, ulcer	10 months	2.5x2.2x1.5	Moderately	T ₂ N ₀ M ₀	HT00625
F10	73	M	Rectum	Adenocarcinoma, ulcer	10 months	4x3	Moderately	T ₂ N ₁ M ₀	HT00626
F11	43	M	Rectum	Adenocarcinoma, ulcer	1 months	3.6x3.5x1.2	Moderately-Well	T ₂ N ₀ M ₀	HT00627

Primary tumor

Tx - Primary tumor cannot be assessed

T0 - No evidence of primary tumor

Tis - Carcinoma in situ; intraepithelial or invasion of lamina propria

T1 - Tumor invades submucosa

T2 - Tumor invades muscularis propria

T3 - Tumor invades through muscularis propria into subserosa or into non-peritonealized pericolic or perirectal tissues.

T4 - Tumor directly invades other organs or structures and/or perforate visceral peritoneum

N - Regional lymph nodes

Nx - Regional lymph nodes cannot be assessed

N0 - No regional lymph node metastasis

N1 - Metastasis in 1 to 3 regional lymph nodes

N2 - Metastasis in 4 or more regional lymph nodes

M - Distant metastasis

Mx - Distant metastasis cannot be assessed

M0 - No distant metastasis

M1 - Distant metastasis

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85