



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

UT243a : Endometrium cancer tissue array with matched adjacent endometrium tissue, including TNM, clinical stage and pathology grade, 12 cases/24 cores, replacing UT243

Microarray Panel	Endometrial endometrioid adenocarcinoma tissue microarray with matched adjacent endometrium tissue, containing 12 cases of endometrium endometrioid adenocarcinoma, duplicate cores per case		
Cores	24		
Cases	12		
Row number	4		
Column number	6		
Core Diameter (mm)	2		
Thickness (µm)	5		
Tissue Array Type	FFPE		
Species	Human		
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.		
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “**” means invalid core; “-” means no applicable or negative in IHC markers.		

	1	2	3	4	5	6
A	Ute	Ute	Ute	Ute	Ute	Ute
B	Ute	Ute	Ute	Ute	Ute	Ute
C	Ute	Ute	Ute	Ute	Ute	Ute
D	Ute	Ute	Ute	Ute	Ute	Ute

Legend: Ute - Uterus

 - Adjacent tissue 1.5 cm away from tumor,
  - Malignant tumor,
  - Malignant tumor (stage I),
  - Malignant tumor (stage IIa),
  - Malignant tumor (stage Ib),
  - Malignant tumor (stage Ic)

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	57	F	Uterus	Endometrioid adenocarcinoma	T1bN0M0	1	Ib	malignant	
A2	2	57	F	Uterus	Cancer adjacent tissue (endometrium tissue s parse)	-	-	-	adjacent	
A3	3	61	F	Uterus	Endometrioid adenocarcinoma	T1N0M0	1	I	malignant	
A4	4	61	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
A5	5	43	F	Uterus	Endometrioid adenocarcinoma	T1N0M0	2	I	malignant	
A6	6	43	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
B1	7	48	F	Uterus	Endometrioid adenocarcinoma	T1bN0M0	1	Ib	malignant	
B2	8	48	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
B3	9	70	F	Uterus	Endometrioid adenocarcinoma	T1bN0M0	1	Ib	malignant	
B4	10	70	F	Uterus	Cancer adjacent tissue (smooth muscle tissu e)	-	-	-	adjacent	
B5	11	70	F	Uterus	Endometrioid adenocarcinoma	T2aN0M0	1	IIa	malignant	
B6	12	70	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	

C1	13	58	F	Uterus	Endometrioid adenocarcinoma	T1bN0M0	1	Ib	malignant	
C2	14	58	F	Uterus	Cancer adjacent tissue (adenocarcinoma s...se)	-	-	-	adjacent	
C3	15	52	F	Uterus	Endometrioid adenocarcinoma	T1cN0M0	2	Ic	malignant	
C4	16	52	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
C5	17	48	F	Uterus	Endometrioid adenocarcinoma	T1cN0M0	2	Ic	malignant	
C6	18	48	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
D1	19	52	F	Uterus	Endometrioid adenocarcinoma	T1bN0M0	2--3	Ib	malignant	
D2	20	52	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
D3	21	61	F	Uterus	Endometrioid adenocarcinoma	T1bN0M0	2--3	Ib	malignant	
D4	22	61	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
D5	23	53	F	Uterus	Endometrioid adenocarcinoma	T1bN0M0	3	Ib	malignant	
D6	24	53	F	Uterus	Cancer adjacent tissue (endometrium tissue)	-	-	-	adjacent	
-	0	42	M	Adrenal gland	Pheochromocytoma (tissue marker)		-		Malignant	

TNM grading:

T - Primary tumor

- Tx - Primary tumor cannot be assessed
- T0 - No evidence of primary tumor
- Tis - Carcinoma in situ; intraepithelial or invasion of lamina propria
- T1 - Tumor invades submucosa
- T2 - Tumor invades muscularis propria
- T3 - Tumor invades through muscularis propria into subserosa or into non-peritonealized pericolic or perirectal tissues.
- T4 - Tumor directly invades other organs or structures and/or perforate visceral peritoneum

N - Regional lymph nodes

- Nx - Regional lymph nodes cannot be assessed
- N0 - No regional lymph node metastasis
- N1 - Metastasis in 1 to 3 regional lymph nodes
- N2 - Metastasis in 4 or more regional lymph nodes

M - Distant metastasis

- Mx - Distant metastasis cannot be assessed
- M0 - No distant metastasis
- M1 - Distant metastasis

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, U.K.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Via Lisano 3
(CP.683)
6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85