



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

UT244a : Chorionic epithelial carcinoma and normal endometrial tissue array, including pathology grade, TNM and clinical stage (AJCC 7.0), 12 cases/24 cores, replacing UT244

Microarray Panel	Chorionic epithelial carcinoma and normal endometrial tissue microarray, containing 10 cases of chorionic epithelial carcinoma, 2 normal endometrial tissue, duplicated cores per cases	
Cores	24	
Cases	12	
Row number	4	
Column number	6	
Core Diameter (mm)	1.5	
Thickness (µm)	5	
Tissue Array Type	FFPE	
Species	Human	
Applications	Routine histology procedures including Immunohistochemistry (IHC) and In Situ Hybridization (ISH), protocols which can be found at our support page.	
Notes	1. TMA slides were sectioned and stored at 4°C and may not be fresh cut, but still suitable for IHC. Please request fresh cut if experiment involves phospho-specific antibodies, RNA studies, FISH or ISH, etc. A minimum of 3 slides per TMA must be purchased to cover the cost of trimming for fresh sectioning. 2. Most TMA slides were not coated with an extra layer of paraffin (tissue cores can be easily seen on the glass). To prevent tissue detachment during antigen retrieval, unbaked slides must be baked for at least 30 to 120 minutes at 60°C. before putting into xylene for de-paraffinization. Baked slides were sent out baked for 2 hours. In the following specsheet, “**” means invalid core; “-” means no applicable or negative in IHC markers.	

	1	2	3	4	5	6
A	Ute	Ute	Ute	Ute	Ute	Ute
B	Ute	Ute	Ute	Ute	Ute	Ute
C	Ute	Ute	Ute	Ute	Ute	Ute
D	Ute	Ute	Ute	Ute	Ute	Ute

Legend: Ute - Uterus
 ● - Malignant tumor (stage IB), ● - Normal tissue

, tissue IDs are available in exported Excel files.

Pos.	No.	Age	Sex	Organ/Anatomic Site	Pathology diagnosis	TNM	Grade	Stage	Type	Image
A1	1	23	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
A2	2	23	F	Uterus	Chorionic epithelial carcinoma (necrosis)	T1N0M0		IB	Malignant	
A3	3	25	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
A4	4	25	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
A5	5	30	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
A6	6	30	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
B1	7	57	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
B2	8	57	F	Uterus	Chorionic epithelial carcinoma (necrosis)	T1N0M0		IB	Malignant	
B3	9	42	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
B4	10	42	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
B5	11	46	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
B6	12	46	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
C1	13	50	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
C2	14	50	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	

C3	15	54	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
C4	16	54	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
C5	17	20	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
C6	18	20	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
D1	19	40	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
D2	20	40	F	Uterus	Chorionic epithelial carcinoma	T1N0M0		IB	Malignant	
D3	21	18	F	Uterus	Endometrial tissue	-		-	Normal	
D4	22	18	F	Uterus	Endometrial tissue	-		-	Normal	
D5	23	42	F	Uterus	Endometrial tissue	-		-	Normal	
D6	24	42	F	Uterus	Endometrial tissue	-		-	Normal	
-	0									

TNM grading:

T - Primary tumor

- Tx - Primary tumor cannot be assessed
- T0 - No evidence of primary tumor
- Tis - Carcinoma in situ; intraepithelial or invasion of lamina propria
- T1 - Tumor invades submucosa
- T2 - Tumor invades muscularis propria
- T3 - Tumor invades through muscularis propria into subserosa or into non-peritonealized pericolic or perirectal tissues.
- T4 - Tumor directly invades other organs or structures and/or perforate visceral peritoneum

N - Regional lymph nodes

- Nx - Regional lymph nodes cannot be assessed
- N0 - No regional lymph node metastasis
- N1 - Metastasis in 1 to 3 regional lymph nodes
- N2 - Metastasis in 4 or more regional lymph nodes

M - Distant metastasis

- Mx - Distant metastasis cannot be assessed
- M0 - No distant metastasis
- M1 - Distant metastasis

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, U.K.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Via Lisano 3
(CP.683)
6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85