



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet

Product Name: Monkey Tissue Array – 22 types of cyno monkey normal tissues covering all major organs


Catalog No.: Z7020003

Lot No.: B404085

Species: Cynomolgus Monkey

Tissue Type: Normal

Tissue Array Diagram:

Cores	Size	Cut	Format	QA/QC	
22	2 mm	4 µm	4X6	H&E, IHC Vimentin/Cytokeratin	

Recommended applications: Whole Body RNA or protein tissue profiling using IHC or ISH for cyno monkey; Antibody characterization

Description: There are 22 tissue types. The tissues were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or APES coated Superfrost slides. They all have a guaranteed six months' shelf-life at 4°C from the date of shipment.

Array position	Sex	Age (yr)	Anatomic site	Histology
A1	M	5	Adrenal gland	Normal
A2	M	3.5	Bladder	Normal
A3	M	10	Cerebellum	Normal
A4	M	10	Cerebral cortex	Normal
A5	M	5.4	GI-Esophagus	Normal
A6	M	5.4	GI-Stomach	Normal
B1	M	5	GI-Small intestine	Normal
B2	M	5	GI-Colon	Normal
B3	M	5.4	Heart	Normal
B4	M	5.4	Kidney cortex	Normal
B5	M	5.4	Kidney medulla	Normal
B6	M	10	Liver	Normal
C1	M	5	Lung	Normal
C2	F	5	Ovary	Normal
C3	M	5	Pancreas	Normal
C4	M	3.5	Skin	Normal
C5	F	6	Spinal cord	Normal
C6	M	5.4	Spleen	Normal
D1	M	5	Striated muscle	Normal
D2	M	5.4	Thymus	Normal
D3	M	5	Thyroid	Normal
D4	F	7	Uterus	Normal

Notes: Bake at 60°C for 30 minutes before use. If antigen retrieving is needed, it is always a good idea to start with a protocol with weak to mild strength.

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: 

AMS BIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
AMS Biotechnology (Europe) Ltd
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
amsbio LLC
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 617 945 5033 or
T: +1 800 987 0985
F: +1 617 945 8218

 **Germany**
AMS Biotechnology (Europe) Ltd
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
AMS Biotechnology (Europe) Ltd
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85