

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

Pig Albumin Antibody



| g / | , and body | | | LA | BORATORIES, INC | LIFE SCIENCES | |
|--|---|---|--|----|-----------------|---------------|--|
| Goat Polyclonal | | Conjugate | HRP | | | | |
| Antigen Affinity | Purified | | | | | | |
| Catalog No. | 100-110P | Uniprot ID | P08835 | | | | |
| Lot No. | 34 | GenelD | 396960 | | | | |
| APPLICATIONS | WB, IHC, ICC, | WB, IHC, ICC, ELISA | | | | | |
| SPECIES REACTIVI | TY Pig | | | | | | |
| AMOUNT | 1 ml | | | | | | |
| CONCENTRATION | 1 mg/ml | 1 mg/ml | | | | | |
| STORAGE/SHELF I | .IFE 2 – 8°C / 1 ye | 2 - 8°C / 1 year from date of receipt | | | | | |
| PHYSICAL STATE | Liquid | Liquid | | | | | |
| BUFFER | Phosphate Bu | Phosphate Buffered Saline (PBS) containing 0.2% BSA and 0.05% Pro-Clean 400 | | | | | |
| ISOTYPE | IgG | IgG | | | | | |
| ORIGIN USA | | | | | | | |
| PRODUCTION PROCEDURES | by affinity chr | Antiserum was solid phase adsorbed to ensure class specificity. The antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to horseradish peroxidase (HRP). | | | | | |
| | | Prior to conjugation, immunoglobulin concentration was determined using Beer's Law where 1mg/mL IgG has an A280 of 1.4. Molar enzyme/antibody protein ratio is 4:1. | | | | | |
| | | By immunoelectrophoresis and ELISA this antibody reacts specifically with pig albumin. This antibody may cross react with albumin from other species. | | | | | |
| APPLICATIONS | | Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. | | | | | |
| | Western Blot Immunohisto Immunocytoo ELISA | chemistry 1:2 chemistry 1:2 | 000 - 1:30,000 200 - 1:500 200 - 1:500 0,000 - 1:50,000 | 1 | | | |
| APPLICATION NO | FES Not all listed a | Not all listed applications have been specifically tested by our laboratory. | | | | | |
| ADDITIONAL INFO https://www.fortislife.com/p/A100-110P Use the link above to view SDS, a current list of citations, and other product specifi | | | | | | information. | |
| | | | | | | | |

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc. Michael Spencer, PhD Date: December 11, 2024

For in vitro laboratory use only. Not for any clinical, therapeutic or diagnostic use in humans or animals. Not for human or animal consumption. This product may not be resold or modified for resale without the prior written approval of Bethyl Laboratories, Inc. The information provided in this data sheet is believed to be correct but does not purport to be all-inclusive and is intended to be used as a guide. Bethyl Laboratories, Inc. shall not be liable or responsible in any way for use of either this information or the material supplied. Disposal of hazardous material may be subject to federal, state or local laws or regulations.