

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

Rabbit IgC Antibody	5 Heav	y and Light Chain Cross-Adsorbed		
Donkey Polycl	onal	Conjugate DyLight® 755		
Antigen Affinity Purified				
Catalog No. A120–208D7				
Lot No.	15			
APPLICATIONS		IHC, ICC, Flow Cyt, IF		
SPECIES REACTIVITY		Rabbit. Minimum reactivity to bovine, chicken, goat, human, mouse, pig and rat		
AMOUNT		1 ml		
CONCENTRATION		0.5 mg/ml		
STORAGE/SHELF LIFE		2 - 8°C / 1 year from date of receipt		
PHYSICAL STATE		Liquid		
BUFFER		Phosphate Buffered Saline (PBS) containing 0.2% BSA and 0.09% Sodium Azide		
FLUOROPHORE/PROTEIN		3.3		
ISOTYPE		lgG		
ORIGIN		USA		
PRODUCTION PROCEDURES		Antiserum was cross adsorbed using bovine, chicken, goat, human, mouse, pig and rat immunosorbents to remove cross reactive antibodies. The antibody to rabbit IgG was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 755.		
		Immunoglobulin concentration was determined using Beer's Law where $1 mg/mL$ IgG has an A280 of 1.4.		
		By immunoelectrophoresis and ELISA this antibody reacts specifically with rabbit IgG and with light chains common to other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to bovine, chicken, goat, human, mouse, pig and rat IgG was detected. This antibody may cross react with IgG from other species.		
APPLICATIONS		Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.		
		Immunohistochemistry 1:50 - 1:500 Immunocytochemistry 1:50 - 1:500 Flow Cytometry 1:50 - 1:200 Immunofluorescence 1:50 - 1:500		
APPLICATION NOTES		Not all listed applications have been specifically tested by our laboratory.		
		DyLight® 755 is excited at 754 (in PBS) and emits at 776 (in PBS).		
ADDITIONAL INFO		DyLight [®] is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. https://www.fortislife.com/p/A120-208D7 Use the link above to view SDS, a current list of citations, and other product specific information.		

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc. Michael Spencer, PhD Date: September 3, 2024