



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Monkey IgG Heavy and Light Chain Antibody

Goat Polyclonal
Antigen Affinity Purified
Catalog No. A140-102A
Lot No. 240119

| | | | | | | | | | |
|-----------------------|---|--------------|-------------------|----------------------|----------------|---------------------|----------------|-------|---|
| APPLICATIONS | WB, IHC, ICC, ELISA | | | | | | | | |
| SPECIES REACTIVITY | Monkey | | | | | | | | |
| AMOUNT | 1 ml | | | | | | | | |
| CONCENTRATION | 1 mg/ml | | | | | | | | |
| STORAGE/SHELF LIFE | 2 – 8°C / 2 years from date of receipt | | | | | | | | |
| PHYSICAL STATE | Liquid | | | | | | | | |
| BUFFER | Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide | | | | | | | | |
| ISOTYPE | IgG | | | | | | | | |
| ORIGIN | USA | | | | | | | | |
| PRODUCTION PROCEDURES | <p>The antibody was isolated by affinity chromatography using antigen coupled to agarose beads.</p> <p>Immunoglobulin concentration was determined using Beer’s Law where 1 mg/mL IgG has an A280 of 1.4.</p> <p>By immunoelectrophoresis and ELISA this antibody reacts specifically with monkey IgG and with light chains common to other monkey immunoglobulins. Species detected are Rhesus, Baboon, Carb Green, Cyno, Chimp and African Green. Reaction to human was also detected. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species.</p> | | | | | | | | |
| APPLICATIONS | <p>Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.</p> <table><tr><td>Western Blot</td><td>1:1000 – 1:30,000</td></tr><tr><td>Immunohistochemistry</td><td>1:200 – 1:2000</td></tr><tr><td>Immunocytochemistry</td><td>1:200 – 1:2000</td></tr><tr><td>ELISA</td><td>1:1000 – 1:30,000; for coating plates 1:100 – 1:500</td></tr></table> | Western Blot | 1:1000 – 1:30,000 | Immunohistochemistry | 1:200 – 1:2000 | Immunocytochemistry | 1:200 – 1:2000 | ELISA | 1:1000 – 1:30,000; for coating plates 1:100 – 1:500 |
| Western Blot | 1:1000 – 1:30,000 | | | | | | | | |
| Immunohistochemistry | 1:200 – 1:2000 | | | | | | | | |
| Immunocytochemistry | 1:200 – 1:2000 | | | | | | | | |
| ELISA | 1:1000 – 1:30,000; for coating plates 1:100 – 1:500 | | | | | | | | |
| APPLICATION NOTES | Not all listed applications have been specifically tested by our laboratory. | | | | | | | | |
| ADDITIONAL INFO | <p>https://www.fortislifesciences.com/p/A140-102A</p> <p>Use the link above to view SDS, a current list of citations, and other product specific information.</p> | | | | | | | | |

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.
Michael Spencer, PhD

Date: January 24, 2024