



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

4-HNE Antibody

Rabbit Polyclonal


Antigen Affinity Purified

Catalog No. A150-121A-0.2MG

Lot No. 0

AMOUNT	0.2 mg
CONCENTRATION	1 mg/ml
STORAGE/SHELF LIFE	2 – 8°C / 1 year from date of receipt
PHYSICAL STATE	Liquid
BUFFER	75 mM Sodium Phosphate, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN ₃ , pH 7.2, 30% glycerol
ISOTYPE	IgG
ORIGIN	USA
APPLICATIONS	Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.
APPLICATION NOTES	Not all listed applications have been specifically tested by our laboratory.
ADDITIONAL INFO	https://www.fortislife.com/p/A150-121A-0.2MG Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.
Michael Spencer, PhD Date: April 5, 2024



Detection of 4-HNE by western blot. *Samples:* 1 μ g BSA conjugated to 4-HNE, Malondialdehyde (MDA), 4-Hydroxy-2-hexenal (4-HHE), Hexanoyl Lysine (H-Lys), and Crotonaldehyde (Croto). *Antibody:* Goat anti-4-HNE antibody (A150-121A-0.2MG) used at 1:1000. *Secondary:* HRP-conjugated donkey anti-goat IgG (A50-101P). *Detection:* Chemiluminescence with an exposure time of 30 seconds.