



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien


T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Detection of 4-HNE by western blot. *Samples:* Whole cell lysate (50 μ g) from Hep-G2 cells untreated (-) or treated with 100 μ M 4-HNE for 1 hr. *Antibody:* Rabbit anti-4-HNE recombinant monoclonal antibody [BLR303M] (A700-303 lot 1) used at 1:1000. *Secondary:* HRP-conjugated goat anti-rabbit IgG (A120-101P). *Detection:* Chemiluminescence with an exposure time of 30 seconds. Lower Panel: Rabbit anti-Actin recombinant monoclonal antibody [BLR057F] (A700-057).



Detection of 4-HNE by western blot. *Samples:* 1 μ g BSA conjugated to 4-HNE, Malondialdehyde (MDA), 4-Hydroxy-2-hexenal (4-HHE), Hexanoyl Lysine (H-Lys), and Crotonaldehyde (Croto). *Antibody:* Rabbit anti-4-HNE recombinant monoclonal antibody [BLR303M] (A700-303 lot 1) used at 1:1000. *Secondary:* HRP-conjugated goat anti-rabbit IgG (A120-101P). *Detection:* Chemiluminescence with an exposure time of 10 seconds.