



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

MDA-BSA

Purified

Catalog No. P150-126-0.5MG

Lot No. 0

AMOUNT	0.5 mg
CONCENTRATION	1 mg/ml
STORAGE/SHELF LIFE	Store at -20°C / 1 year from date of receipt
PHYSICAL STATE	Liquid
BUFFER	10 mM PBS, 0.15 M NaCl, 0.5 mM EDTA, 0.01% NaN ₃ , 50% glycerol
ISOTYPE	IgG
ORIGIN	USA
APPLICATIONS	Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.

ADDITIONAL INFO <https://www.fortislife.com/p/P150-126-0.5MG>
Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.
Michael Spencer, PhD Date: April 5, 2024