



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



---

## PBS Buffer Tablets

---

[nordicmubio.com/product/pbs-buffer-tablets](http://nordicmubio.com/product/pbs-buffer-tablets)

Catalogue number: **PU000**

Product Type	Buffers and Reagents
Units	10 L
Application	ELISA Immunohistochemistry Western Blotting

### Background

10 buffer tablets in 1000 ml distilled water give a 10 mM PBS solution pH 7.3. Phosphate buffer tablets for isotonic PBS phosphate buffered saline for universal application in IHC, WB, ELISA

### Product

100 tablets

*Purification Method:* 100 tablets

### Applications

IHC, ELISA , WB

*Working Concentration:* (100 tablets) 1 tablet per 100 ml distilled water

### Storage

2-25°C

### Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. It may contain hazardous ingredients. Please refer to the Safety Data Sheets (SDS) for additional information and proper handling procedures. Dispose product remainders according to local regulations. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.

### References

