



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Product Information

<b>Description:</b>	This liquid Assay Medium is optimized for use in cellular assays in select BPS Bioscience cell lines.
<b>Formulation:</b>	<ul style="list-style-type: none"><li>• <b>Assay Medium 5A (BPS Bioscience #79620):</b> Opti-MEM Reduced Serum medium (Thermo Fisher, #31985070), 0.5% calf serum (Cytiva, #SH30072.04), 1% non-essential amino acids (Corning, #25-025-CI), 1 mM Na pyruvate (Corning, #25-000-CI), 10 mM HEPES (Cytiva SH30237.01), 1% Penicillin/Streptomycin (Thermo Fisher, #15140163).</li></ul>
<b>Quality:</b>	Sterile filtered (0.2 $\mu$ m) and Mycoplasma-free. Each lot is tested for growth of recombinant cells. All components are low endotoxin. FBS is ISIA-certified and sourced from the United States of America.
<b>Storage:</b>	Upon receipt, the medium should be stored in the dark at 4°C.
<b>Stability:</b>	Do not exceed 37°C. When stored in the dark at 4°C, the product is stable for six months after first thaw. Avoid repeated freeze/thaw cycles.
<b>Culture Protocol:</b>	For detailed cellular assay protocol, see cell line data sheet.