



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.202.1401  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## **Thaw Medium 11**

### **Catalog #: 79976**

#### **Description**

This liquid medium is optimized to culture select BPS Bioscience cell lines.

#### **Formulation**

DMEM medium (Hyclone #SH30024.01) supplemented with 10% heat-inactivated FBS (Gibco #26140-079), 1% GlutaMAX (Gibco #35050-061), 1% Penicillin/Streptomycin (Hyclone #SV30010.01).

#### **Quality**

Sterile filtered (0.2  $\mu$ m) and Mycoplasma-free. Each lot is tested for growth of recombinant cells. All components are low endotoxin. FBS is ISIA-certified and sourced from the United States of America.

#### **Storage and Stability**

Upon receipt, the medium should be stored in the dark at 4°C. Avoid repeated freeze/thaw cycles. Do not exceed 37°C. When stored in the dark at 4°C, the product is stable for six months after first thaw.

#### **Culture Protocol**

For detailed thawing and cell culture protocol, see cell line data sheet.

#### **Optimized for Use With**

<b>Product</b>	<b>Cat #</b>	<b>Size</b>
NF- $\kappa$ B Reporter (Luc) – Raw 264.7 Cell line	#79978	2 vials

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)