



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Product Information

<b>Description:</b>	The NP-40 Cell Lysis Buffer is intended for the lysis of mammalian adherent and suspension-cultured cells. This buffer is suitable for extraction of cytoplasmic proteins in mild, non-denaturing conditions.
<b>Formulation:</b>	50 mM Tris HCl, pH 7.4, 150 mM NaCl, 1 mM EDTA, 1 mM EGTA, and 1% IGEPAL CA-630
<b>Quality:</b>	Sterile filtered (0.2 $\mu$ m)
<b>Storage/Stability:</b>	Upon receipt, store at 4°C. Stable for 6 months from date of receipt, when stored as directed.
<b>Instructions for Use:</b>	This is a ready-to-use buffer. Mix well before use. This buffer does not contain protease or phosphatase inhibitors. If desired, add appropriate inhibitors (such as protease inhibitors, phosphatase inhibitors, ADP-ribosylation inhibitors etc.) immediately before use to prevent protein proteolysis or to maintain phosphorylation and/or ribosylation status of proteins.
<b>Shipping Temperature:</b>	4°C (gel packs)