

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Cyclosporin A

Item No. 12088

CAS Registry No.: 59865-13-3 NSC 290193 Synonym: MF: $C_{62}H_{111}N_{11}O_{12}$

FW: **Purity:** ≥99% Supplied as: A solid -20°C Storage: Stability: ≥2 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Cyclosporin A is supplied as a solid. A stock solution may be made by dissolving the cyclosporin A in the solvent of choice, which should be purged with an inert gas. Cyclosporin A is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF). The solubility of cyclosporin A in these solvents is approximately 14, 3, and 20 mg/ml, respectively.

Cyclosporin A is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, cyclosporin A should first be dissolved in DMF and then diluted with the aqueous buffer of choice. Cyclosporin A has a solubility of approximately 0.25 mg/ml in a 1:3 solution of DMF:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

Cyclosporin A is an immunosuppressant drug that binds cyclophilin D, inhibiting the phosphatase activity of calcineurin in T cells ($IC_{50} = 5 \text{ nM}$). 1-3 It blocks the transcription of genes involved in T cell activation by interfering with the calcineurin-driven nuclear import of the transcription factor NFAT.^{1,4,5} Because of its powerful immunosuppressant activity, cyclosporin A has long been used for the prevention of transplant rejection.6

References

- 1. Jurgensen, K.A., Koefoed-Nielsen, P.B., and Karamperis, N. Calcineurin phosphatase activity and immunosuppression. A review on the role of calcineurin phosphatase activity and the immunosuppressive effect of cyclosporin A and tacrolimus. Scand. J. Immunol. 57(2), 93-98 (2003).
- 2. Mrowietz, U. Macrolide immunosuppressants. Eur. J. Dermatol. 9(5), 346-351 (1999).
- Fruman, D.A., Klee, C.B., Bierer, B.E., et al. Calcineurin phosphatase activity in T lymphocytes is inhibited by FK 506 and cyclosporin A. Proc. Natl. Acad. Sci. USA 89(9), 3686-3690 (1992).
- 4. Roehrl, M.H.A., Kang, S., Aramburu, J., et al. Selective inhibition of calcineurin-NFAT signaling by blocking protein-protein interaction with small organic molecules. Proc. Natl. Acad. Sci. USA 101(20), 7554-7559
- 5. McCaffrey, P.G., Perrino, B.A., Soderling, T.R., et al. NF-ATp, a T lymphocyte DNA-binding protein that is a target for calcineurin and immunosuppressive drugs. J. Biol. Chem. 268(5), 3747-3752 (1993).
- Laupacis, A., Keown, P.A., Ulan, R.A., et al. Cyclosporin A: A powerful immunosuppressant. Can. Med. Assoc. J. 126(9), 1041-1046 (1982).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 01/03/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM