

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# PRODUCT INFORMATION



## **BI-D1870**

Item No. 15264

CAS Registry No.: 501437-28-1

Formal Name: 2-[(3,5-difluoro-4-hydroxyphenyl)

amino]-7,8-dihydro-5,7-dimethyl-8-(3-

methylbutyl)-6(5H)-pteridinone

MF:  $C_{19}H_{23}F_2N_5O_2$ 

FW: 391.4 ≥95% **Purity:** 

UV/Vis.:  $\lambda_{\text{max}}$ : 211, 279 nmm Supplied as: A crystalline solid

Storage:

As supplied, 2 years from the QC date provided on the Certificate of Analysis, when Stability:

stored properly

#### **Laboratory Procedures**

BI-D1870 is supplied as a crystalline solid. A stock solution may be made by dissolving the BI-D1870 in the solvent of choice. BI-D1870 is soluble in organic solvents such as DMSO and dimethyl formamide, which should be purged with an inert gas. The solubility of BI-D1870 in these solvents is approximately 15 and 10 mg/ml, respectively.

BI-D1870 is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, BI-D1870 should first be dissolved in DMSO and then diluted with the aqueous buffer of choice. BI-D1870 has a solubility of approximately 0.3 mg/ml in a 1:2 solution of DMSO:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

#### Description

The 90 kDa ribosomal S6 kinases (RSKs) are a group of serine/threonine kinases involved in diverse cellular processes, including growth, survival, and motility. BI-D1870 is a cell permeable, ATP-competitive inhibitor of the four vertebrate isoforms of RSK, RSK1-4 (IC<sub>50</sub>s = 31, 24, 18, and 15 nM, respectively).<sup>2</sup> At 100 nM, it also significantly inhibits polo-like kinase 1, Aurora B, maternal embryonic leucine zipper kinase, and mammalian STE20-like kinase 2.2,3 BI-D1870 acts at the N-terminal kinase domain, but not the C-terminal domain, of RSK.<sup>2,3</sup>

#### References

- 1. Romeo, Y., Zhang, X., and Roux, P.P. Regulation and function of the RSK family of protein kinases. Biochem. J. 441(2), 553-569 (2012).
- 2. Sapkota, G.P., Cummings, L., Newell, F.S., et al. BI-D1870 is a specific inhibitor of the p90 RSK (ribosomal S6 kinase) isoforms in vitro and in vivo. Biochem. J. 401(1), 29-38 (2007).
- Bain, J., Plater, L., Elliot, M., et al. The selectivity of protein kinase inhibitors: A further update. Biochem. J. 408, 297-315 (2007).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

### WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 03/23/2016

#### **CAYMAN CHEMICAL**

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM