

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



6-diazo-5-oxo-L-nor-Leucine

Item No. 17580

CAS Registry No.: 157-03-9

Formal Name: 6-diazo-5-oxo-L-norleucine

Synonyms: **DON, NSC 7365** MF: $C_6H_9N_3O_3$ FW: 171.2 **Purity:** ≥95%

 λ_{max} : 274 nm A crystalline solid UV/Vis.: Supplied as:

Storage: -20°C Stability: ≥2 vears

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

6-diazo-5-oxo-L-nor-Leucine (DON) is supplied as a crystalline solid. A stock solution may be made by dissolving the DON in the solvent of choice, which should be purged with an inert gas. DON is soluble in the organic solvent DMSO at a concentration of approximately 1 mg/ml.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of DON can be prepared by directly dissolving the crystalline solid in aqueous buffers. The solubility of DON in PBS, pH 7.2, is approximately 10 mg/ml. We do not recommend storing the aqueous solution for more than one day.

Description

DON is a glutamine analog that inhibits glutaminases ($K_i = 6 \mu M$) and other glutamine-utilizing enzymes, including cytidine triphosphate synthase (CTPS), which uses glutamine in the synthesis of CTP.1-3 DON is commonly used in cells or tissues but its use in whole animals is limited by its teratogenicity.³⁻⁵

References

- 1. Willis, R.C. and Seegmiller, J.E. The inhibition by 6-diazo-5-oxo-I-norleucine of glutamine catabolism of the cultured human lymphoblast. J. Cell Physiol. 93, 375-382 (1977).
- Thangavelu, K., Chong, Q.Y., Low, B.C., et al. Structural basis for the active site inhibition mechanism of human kidney-type glutaminase (KGA). Sci. Rep. 4, 1-7 (2014).
- Thomas, A.G., Rojas, C., Tanega, C., et al. Kinetic characterization of ebselen, chelerythrine and apomorphine as glutaminase inhibitors. Biochem. Biophys. Res. Commun. 438(2), 243-248 (2013).
- Botman, D., Tigchelaar, W., and Van Noorden, C.J.F. Determination of phosphate-activated glutaminase activity and its kinetics in mouse tissues using metabolic mapping (quantitative enzyme histochemistry). J. Histochem. Cytochem. 62(11), 813-826 (2014).
- 5. Greene, R.M. and Kochhar, D.M. Limb development in mouse embryos: Protection against teratogenic effects of 6-diazo-5-ox-L-norleucine (DON) in vivo and in vitro. J. Embryol. Exp. Morphol. 33(2), 355-370 (1975).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 06/09/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM