

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# PRODUCT INFORMATION



## Flibanserin

Item No. 19203

CAS Registry No.: 167933-07-5

Formal Name: 1,3-dihydro-1-[2-[4-[3-(trifluoromethyl)phenyl]-

1-piperazinyl]ethyl]-2H-benzimidazol-2-one

Synonyms: BIMT 17, BIMT 17BS

MF:  $C_{20}H_{21}F_3N_4O$ 

FW: 390.4 **Purity:** ≥98%

Supplied as: A crystalline solid

Storage: -20°C

Stability: As supplied, 2 years from the OC date provided on the Certificate of Analysis, when

stored properly

## **Laboratory Procedures**

Flibanserin is supplied as a crystalline solid. A stock solution may be made by dissolving the flibanserin in the solvent of choice. Flibanserin is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF), which should be purged with an inert gas. The solubility of flibanserin in ethanol is approximately 15 mg/ml, and 30 mg/ml in DMSO and DMF.

Flibanserin is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, flibanserin should first be dissolved in DMF and then diluted with the aqueous buffer of choice. Flibanserin has a solubility of approximately 0.2 mg/ml in a 1:4 solution of DMF:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

#### Description

Flibanserin is a full agonist of the serotonin 5-HT $_{1A}$  receptor and an antagonist of 5-HT $_{2A}$  (K $_{i}$ s = 1 and 49 nM, respectively). It also binds to dopamine D $_{4}$  receptors with K $_{i}$  values ranging from 4-24 nM, but demonstrates no affinity for the other 5-HT subtypes or other neurotransmitter receptors. In vitro, flibanserin has been shown to reduce forskolin-stimulated cAMP formation in cells and rat tissues and to antagonize the accumulation of phosphatidyl inositol turnover induced by 5-HT in the mouse cortex. At 10 mg/kg, flibanserin can reduce serotonin in the prefrontal cortex and dorsal raphe of conscious rats while increasing extracellular noradrenaline and dopamine.1 Multifunctional serotonergic ligands like flibanserin that can enhance downstream release of dopamine and norepinephrine while concomitantly reducing serotonin release in the brain circuits that mediate symptoms of reduced sexual interest and desire are being studied clinically for therapeutic potential to improve sexual functioning.<sup>2</sup>

#### References

- 1. Invernizzi, R.W., Sacchetti, G., Parini, S., et al. Flibanserin, a potential antidepressant drug, lowers 5-HT and raises dopamine and noradrenaline in the rat prefrontal cortex dialysate: Role of 5-HT<sub>1Δ</sub> receptors. Br. J. Pharmacol. 139(7), 1281-1288 (2003).
- 2. Stahl, S.M. Mechanism of action of flibanserin, a multifunctional serotonin agonist and antagonist (MSAA), in hypoactive sexual desire disorder. CNS Spectr. 20(1), (2015).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

al should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 05/02/2016

## **CAYMAN CHEMICAL**

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM