



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

PRODUCT INFORMATION



Calcein Red™ AM

Item No. 20632

Synonym: Calcein Red™ Acetoxymethyl ester
FW: 1,015.8
Purity: ≥90%
Ex./Em. Max: 560/574 nm
Supplied as: A solid
Storage: -20°C
Stability: ≥2 years
Special Conditions: Light sensitive

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Calcein Red™ AM is supplied as a solid. A stock solution may be made by dissolving the Calcein Red™ AM in the solvent of choice. Calcein Red™ AM is soluble in the organic solvent DMSO, which should be purged with an inert gas.

Description

Calcein Red™ AM is a fluorogenic dye that is used to assess cell viability. As the acetoxymethyl ester, Calcein Red™ AM is non-fluorescent but cell-permeant and is hydrolyzed to the fluorescent probe Calcein Red™ (Item No. 20633) by intracellular esterases. Calcein Red™ is retained by living cells and displays excitation/emission maxima of 560/574 nm, respectively. Calcein Red™ AM is a red-emitting variant of the green fluorophore calcein AM (Item No. 14948), which has been used to evaluate cell viability and cytotoxicity.^{1,2}

References

1. Bratosin, D., Mitrofan, L., Palii, C., *et al.* Novel fluorescence assay using calcein-AM for the determination of human erythrocyte viability and aging. *Cytometry* **66**(1), 78-84 (2005).
2. Lévesque, A., Paquet, A., and Pagé, M. Measurement of tumor necrosis factor activity by flow cytometry. *Cytometry* **20**(2), 181-184 (1995).

WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 03/23/2018

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897
[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM