

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# PRODUCT INFORMATION



## **Thiothixene**

Item No. 23649

CAS Registry No.: 5591-45-7

N,N-dimethyl-9-[3-(4-methyl-1-piperazinyl) Formal Name:

propylidene]-9H-thioxanthene-2-sulfonamide

Synonyms: CP 12,252-1, NSC 108165, P 4657B,

cis-Thiothixene

MF:  $C_{23}H_{29}N_3O_2S_2$ 

FW: 443.6 **Purity:** 

UV/Vis.:  $\lambda_{\text{max}}$ : 230, 307 nm Supplied as: A crystalline solid

-20°C Storage: Stability: ≥4 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



Thiothixene is supplied as a crystalline solid. A stock solution may be made by dissolving the thiothixene in the solvent of choice, which should be purged with an inert gas. Thiothixene is soluble in organic solvents such as DMSO and dimethyl formamide. The solubility of thiothixene in these solvents is approximately 0.2 and 0.5 mg/ml, respectively.

### Description

Thiothixene is a typical antipsychotic.  $^1$  It selectively binds to dopamine  $D_2$  over  $D_1$ ,  $D_3$ , and  $D_4$  receptors (Kis = 0.417, 338, 186.2, and 363.1 nM, respectively). Thiothixene also binds to various serotonin (5-HT), histamine  $H_1$ ,  $\alpha_1$ - and  $\alpha_2$ -adrenergic, muscarinic acetylcholine, and sigma receptors (K.s = 15-5,754 nM) as well as the dopamine, norepinephrine, and serotonin transporters (K<sub>i</sub>s = 3.16-30 μM). In vivo, thiothixene reduces spontaneous and amphetamine-induced locomotor activity in rats.<sup>2</sup> It enhances latent inhibition, as measured by a decreased lick latency in response to light and foot shock stimuli, which is a measure of selective attention in rats.<sup>3</sup> Thiothixene also increases competitive behavior in submissive mice, indicating antidepressant-like behavior.<sup>4</sup> Formulations containing thiothixene have been used in the treatment of schizophrenia and bipolar mania.

#### References

- 1. Silvestre, J.S. and Prous, J. Research on adverse drug events. I. Muscarinic M<sub>3</sub> receptor binding affinity could predict the risk of antipsychotics to induce type 2 diabetes. Methods Find. Exp. Clin. Pharmacol. 27(5), 289-304 (2005).
- 2. Schaefer, G.J. and Michael, R.P. Drug interactions on spontaneous locomotor activity in rats. Neuroleptics and amphetamine-induced hyperactivity. Neuropharmacology 23(8), 909-914 (1984).
- Dunn, L.A., Atwater, G.E., and Kilts, C.D. Effects of antipsychotic drugs on latent inhibition: Sensitivity and specificity of an animal behavioral model of clinical drug action. Psychopharmacology (Berl). 112(2-3), 315-323 (1993).
- 4. Malatynska, E., Rapp, R., Harrawood, D., et al. Submissive behavior in mice as a test for antidepressant drug activity. Pharmacol. Biochem. Behav. 82(2), 306-313 (2005).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 11/02/2022

#### **CAYMAN CHEMICAL**

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM