

## Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

#### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

# **PRODUCT** INFORMATION



**D-Desthiobiotin** 

Item No. 30418

CAS Registry No.:	533-48-2	
Formal Name:	(4R,5S)-5-methyl-2-oxo-4-	Н
	imidazolidinehexanoic acid	
MF:	C <sub>10</sub> H <sub>18</sub> N <sub>2</sub> O <sub>3</sub>	N
FW:	214.3	
Purity:	≥95%	N OH
Supplied as:	A crystalline solid	/ с
Storage:	-20°C	
Stability:	≥4 years	

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

#### Laboratory Procedures

D-Desthiobiotin is supplied as a crystalline solid. A stock solution may be made by dissolving the D-desthiobiotin in the solvent of choice, which should be purged with an inert gas. D-Desthiobiotin is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF). The solubility of D-desthiobiotin in ethanol is approximately 0.5 mg/ml and approximately 20 mg/ml in DMSO and DMF.

D-Desthiobiotin is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, D-desthiobiotin should first be dissolved in DMSO and then diluted with the aqueous buffer of choice. D-Desthiobiotin has a solubility of approximately 0.2 mg/ml in a 1:2 solution of DMSO:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

#### Description

D-Desthiobiotin is a protein cross-linking agent and a precursor in the synthesis of biotin (Item No. 22582).<sup>1</sup> It has been used in the preparation of agarose matrices for affinity-based isolation of streptavidin-fluorophore conjugates.

#### Reference

1. Hirsch, J.D., Eslamizar, L., Filanoski, B.J., et al. Easily reversible desthiobiotin binding to streptavidin, avidin, and other biotin-binding proteins: Uses for protein labeling, detection, and isolation. Anal. Biochem. 308(2), 343-357 (2002).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

#### SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 12/21/2022

#### CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM